

Product datasheet for **SC306184**

ZNF570 (NM_144694) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF570 (NM_144694) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF570
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_144694, the custom clone sequence may differ by one or more nucleotides <pre> ATGGCTGTTGGGCTTCTTAAAGCCATGTACCAGGAGTTGGTGACCTTCAGAGATGTGGCT GTAGACTTCTCCCAAGAGGAATGGGATTGTCTGGATTCTTCTCAAAGACATCTGTACAGT AATGTGATGCTAGAGAACTACAGGATCTTGGTATCACTGGGACTTTGCTTTTCCAAACCA AGTGTGATATTATTGTTGGAACAAGGAAAAGCACCCCTGGATGGTGAAGAGAGAGCTGACA AAAGGCTTGCTCAGGCTGGGAGCCTATATGTGAGACTGAAGAATTAACCCCAAAGCAG GATTTTTATGAAGAACATCAATCCCAGAAGATAATAGAAACACTTACAAGCTATAACCTT GAATACTCCAGTTTGAGAGAAGAGTGGAATGTGAGGGCTATTTTGAAAGGCAACCAGGT AATCAGAAAGCGTGTTCAGGAAGAGATAATCACTCATGAAGAACCCTTTTGTATGAG AGAGAACAAGAATATAAATCTTGGGAAGTTTTCATCAGAACCCACTGCTTGTACACAA AAGATAATCCCCAAAGAGGAGAAAGTACATAAACATGACACACAAAAGAGAAGCTTTAA AAAAATTAATGGCTATTAAGCCCAAGAGTGTCTGTGCAGAGAAGAACTTTGAAATGT AATGACTGTGAAAAAGTCTTCAGCCAGAGTTCATCCCTTACTCTTCATCAAAGAATTCAT ACTGGAGAGAAACCCTATAAATGTATAGAGTGTGGAAGCCCTTCAGCCAGAGATCAAAT CTTGTTCAACATCAGAGGATTCATACTGGAGAAAAACCCTATGAATGTAAGGAATGTAGG AAAGCCTTCAGTCAGAATGCACACCTAGTTCAACATCTGCGAGTTCATACTGGAGAAAA CCTTACGAATGTAAGGTATGTCGAAAAGCCTTCAGCCAGTTTGCTACCTTGCTCAACAT CAGAGAGTTCACACGGGAGAGAAACCCTATGAATGTATCGAATGTGGGAAAGCATTTAGC AACAGATCATCCATTGCTCAACACCAGAGAGTTCATACAGGAGAGAAACCCTATGAATGT AATGTCTGTGGGAAAGCATTTAGCCTTCGTGCATACCTTACTGTACATCAGAGAAACAT ACTGGAGAGAGACCCTATGAATGTAAGGAATGTGGGAAGGCCCTTAGCCAGAATTCACAC CTTGCTCAACATCAGAGAATTCATACTGGAGAAAAACCTTATAAGTGTGAGGAATGTAGG AAAGCATTGAGCCAGATTGCCTACCTTGCTCAGCATCAAAGAGTTCATACTGGAGAGAAA CCCTATGAATGTATTGAATGTGGGAAGGCTTTTAGCAATGACTCGTCCCTTACTCAACAT CAGCGAGTTCATACTGGAGAGAAACCTTATGAATGTACTGTTTGTGGAAGGCTTTTAGT TACTGTGGATCCCTTGCCCAACATCAGAGAATTCATACTGGAGAGAGACCCTATGAATGT AAGGAATGCAAAAAAACCTTCAGGCAGCATGCACACCTTGCTCATCACCAGAGAATTCAC ATTGGGAGTCACTGTCACCAACCAACCCAGTCAATCACCAAGTCCTATAG </pre>
Restriction Sites:	Please inquire
ACCN:	NM_144694



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_144694.1, NP_653295.1</u>
RefSeq Size:	2247 bp
RefSeq ORF:	1611 bp
Locus ID:	148268
UniProt ID:	<u>Q96NI8</u>
Cytogenetics:	19q13.12
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) differs in the 5' UTR, lacks a portion of the 5' coding region and initiates translation at a downstream start codon, compared to variant 1. It encodes isoform 2, which has a shorter N-terminus, compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.