

## Product datasheet for **SC306202**

### CCDC138 (NM\_144978) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	CCDC138 (NM_144978) Human Untagged Clone
Tag:	Tag Free
Symbol:	CCDC138
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC306202 representing NM\_144978.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGAGCCGAGGGTCGTCAAGCCACCGGGCAGGATTTAGTAGTGGAGAGTCTCAAAGCCGCTACGGA
CTCGGGGGCAGCTGCCCGACGAGTATGATTTTTCAAATTTTTATCAGTCTAAGTATAAGAGAAGAACT
CTAACCTCCCCAGGTGATTTGGATATCTACTCTGGAGATAAAGTTGGTTCATCGTTAAAAATTCTGAT
GAAAGCAAGCATTGTAGAACACCATTGGGCAGCTTATTCAAGCACGTAATGTGAATTGCCTAGATGAT
GAACTGGATTCTTTCCATGATTTGAAGAAACAGGAAACAGAAGAAGGTTAATTGAAAATGATTATAGA
GTTAGTACCTCGAAAATAACCAAGCAGTCTTTAAAGAAATAGAAAAAGTTGCCTTGCCAACTAATACG
ACCTCATCGAGACCTCGGACTGAGTGTGTAGTGTGCAGGTGACTCTCCTTTGAAACCTGTGAGTGT
CCAAAACTAAAGCATCAGACAAGCGGAGTTTACTTCCACATCAGATCAGTCAGATATATGACGAATTA
TTTCAGATACATCTGAAATTGCAGTGTGAAACTGCAGCACAACAGAAATTTGCTGAAGAAGTTCAAAG
CGAGAAGCTTTTTACTTGAAAGAGAACAAGTCTTTTTCAGACATGAAAATGCCTTGAGTAAAATTTAA
GGTGTGAAGAAGAGTTCTTACAAGATTTCAAATTTATAAAGAGCAGCATGATGCAGAAGTTGAACAC
TTAACCGAAGTTCTTAAGGAAAAGAATAAAGAAACCAAGAGACTGAGGTCTCTTTTGATGCATTGAAA
GAATTGAATGATACCTTAAAAAACAGTTAAATGAAGCAAGTGAAGAAAACAGGAAGATAGACATTCAG
GCTAAAAGAGTTCAAGCTCGTTTAGATAATTTACAGAGGAAGTACGAGTTTATGACAATACAGAGATTG
AAAGGAAGTTCCCATGCTGTTTATGAAAAGTAAAACAAGAAAAGCACCAGTTTCAAAGACT
TACAAGGTACCACTTAAATGGGCAAGTTTATGAACCTTTAACTGTCTTATGGACTGGATTTCCGATCAT
CATCTTAGCAAAGTGAACATGAAGAATCTGGAATGGATGGTAAAAACCACAAGTCAAATTTGCTTCC
CAGAGAAATGATATTCAGGAGAAGTGTGTAAGCTTTTGCCTCTAATGACAGAGCAGCTACAGTGGATG
CCATTTGTGAATATCAAAGTTCAGGAGCCTTTTGTAAAATTTATATATTGGTCCCTAAGGCAGCTAGAT
GCTGGAGCACAGCACTCGACTATGACATCAACATTGAGGAGATTGGGTGAAGACATTTTTAAAGGAGTG
GTAACAAAGGAATTCAGGATAATTCTCCACAGCATTCTGTGGAGAATAAACCAAGACAGCTGCTTTC
TTAAGAGCTCCAATTTGCCATTGAGATTTTTATCAACCTTAAATGTTCTCAAACAGTCACTCAAGCT
GATTACCTGGCTCAGGCATTTGATTCTCTTTGTTGGACTTGAAGACAGAAGAAGGAAAAACCTTGTTT
TTGGAGTATCAGGCTGTTCCAGTAAATTAAGTCATCTAAGAATATCCAGTAAAGGACTCCTGTCTAAT
GTTATTGATAGTTTGTCCAGATGACGGTGAATCTAAATCCTTGACGCTTTTCCGGAAGCCTGTAGC
AACTCTTTATTTTTCGTACTTGTCTGTGCTGCTTCGAGCCCTAAGCTTGATCTTCAAATACTAGAA
AAACTCAGTATTATTTTACAGAACTTTCCAAAATCAAGAGTAATAAGAAGCTCTTTGAACTTTTTACG
ATTCATCTGATGTTCAAGAAATACAAAGGACAACAAACCCAGAGCATGCATTTCTCTGTATTAATCTA
AATTCAACTCTGTTCAATCTGGGTTTAAACAAAATGTAACCTCCCTGGTCTCCAGTGAAGCCCTTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_144978  
**Insert Size:** 1998 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_144978.2](#)

**RefSeq Size:** 2420 bp

**RefSeq ORF:** 1998 bp

**Locus ID:** 165055

**UniProt ID:** [Q96M89](#)

**Cytogenetics:** 2q13

**MW:** 76.2 kDa