

Product datasheet for **SC306211**

DRC1 (NM_145038) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DRC1 (NM_145038) Human Untagged Clone
Tag:	Tag Free
Symbol:	DRC1
Synonyms:	C2orf39; CCDC164; CILD21
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC306211 representing NM_145038.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGAATCCGCCGGGGTCCCTAGAGGCCCTGGACCCGAACGTGGACGAGCACTTGTCCACCCAGATTCTC
GCGCCCTCGGTCCACTCCGACAACCTCAGGAGCGCATCCAGGCCCGCGCCTCCGCATCGCTGCGCGC
TTAGAAGCCCGGAGCGGGAAGCACTTGGAGAATATTTAGATGGGAAGAAGGAGAGTGAGGAGGATCAA
AGCAAGAGCTACAACAGAAAGAAGAAAGCCGATTGAAATTGGCTAAACTTCTGTCTGTGGCACCGAG
TTGGTGACAAATATCCAGGTGGCTATAGATATCAGAGAGATTACAGGAGAGTCGAAGAAGAGGAGATA
AAGCGTCAAAGGATTGAGAAGCTGGAGAATGAAGTTAAGACCAGCCAGGACAAAATTCGATGAAATCACC
TCAAAGTGGGAAGAGGGCAAGCAGAAGAGAATCCCAAGAGCTGTGGGAAATGCTCAATACCCAACAG
CTGCACGTGCTGGACTTTAGAAGATAAGAATAAACTCATCAGCGAGTTACAGCAGGAGTTAAAAACA
AAGGATGACCAGTATGTGAAGGATTTGAAGAAACAGTCAGATGACATCTGCCTGCTTCTGGAGCGGATG
GAAGAACAGGTGAAGAATGTGATGAAAACCTTTCTGTGAGGAGCTCTATAACATTGAGAAAGCATTGAG
GTGGAACGCCAAGAGCTACTGGCCAGTAATAAGAAGAAATGGGAGCAAGCCCTTACAGGCTATAATGCC
AAAGAGCTGGAGTATCTTAAACAACCGCATGAAGAAAGTAGAGGACTATGAGAAGCAGCTGAACAGGCAG
AGGATCTGGGATTGCGAAGAATACAACATGATCAAGATCAAGCTGGAGCAGGATGTGCAGATTCTTGAG
CAGCAGCTTCAGCAGAGGAAAGCAATTTATCAGCTAAACCAAGAGAAATTAGAGTACAACCTGCAGGTG
CTGAAGAAGAGAGATGAAGAAAGCACAGTAATTAATCCCAGCAGAAGAGGAAGTCAATCGCCTGCAT
GATATACTTAAACAATCTGAGATCAAAATATGCCAAGCAATAAAGCAGTTTCAGGAGGAGAACAGTCT
CTAACCTCGGACTACAACGCTTGTGATGCAATTCAGGAGCTACAGAAAGCCATGAGGCATTTTGCT
CTTATTGATGATGAGAAGTTTTGGGAGATTTGGCTGATGAATGAAGAGGAGGCGAAGGACCTAATAGCC
AGAGCCTTTGATGTGGACAGGATCATCCACCCATCATCTGGGGCTTCCCTGGGCAGCACCTGATTTTC
TGTTCTGAACAATGTTGGCCTATTTCTCAGCAGCCCCAGAAGTCCGCCACACAGATAGTAGAAGAA
ATGCTTATGCGCTCAGAAGAGGAGGAGGCAGAAGAGGCCCGCGGAACAGAGTCTACCTGGATTTG
CCGAAGCAAAATTTCTGAAAAAATACCAAGAGGATCCTGATGCTCCTGTGTGACGAGTCGGGGTTCTC
ATAGAGAGCAAGCTGCTGAGCCTTCTCCTGCCCTGGAGCAGAATGAATGCTATCTGCTGAGGCTGGAT
GCCATCTTCTCCGCCCTTGAATTGAAAGTGAGGATGACTTGTATAAACTGGTGAACCTTCTCCTAAA
TATCGAGCTCACGTTTATCTTCCAGCCTCCAGATCAAGCCCTGCAGTCAGGCGAGCATGGAGAAGGCG
AGCATGGAGGAGACAAGCAGAGGTGCAATTGGAGCTGGCAGAGCAGACGGAGATGGAGGGAGAAAAG
GAAGAAAGCCTGGTGAAGGGGAGAAGGAGGAAGAGGAGGAGACCCACCCTCCCCTGGGTATCCAC
CCCAATGATGTCCTCAAGATTCTGGAGGCCCTCGTCATGGGTCTGAAGAAGCCTAGGGACTCGCGGGCC
CCGCTGAGGGTACAGAAGAATGTGCGTGACAACCTCAAGGACTCGGAGTACTGGCAGGCCCTGACCACA
GTGATCCCTTCTCCAAGCAGAACCTCTGGGATGCCCTCTACACAGCCTTGGAGAAGTACCACCTTGTG
CTGACCCAGAGGGCCAAGCTGCTGCTGAAAAACAGTTCTCTGGAGCAGCAGAACACAGAGCTGCAGGCG
CTACTGCAGCAGTATCTGAACCTCAAGATCAACTCTGAACTGCAAGTTCTCCACTCAGGTGTTGCGG
GTACCCACAAAATGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_145038

Insert Size: 2223 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_145038.4
RefSeq Size:	2519 bp
RefSeq ORF:	2223 bp
Locus ID:	92749
UniProt ID:	Q96MC2
Cytogenetics:	2p23.3
MW:	87.1 kDa
Gene Summary:	This gene encodes a central component of the nexin-dynein complex (N-DRC), which regulates the assembly of ciliary dynein. Mutations in this gene can cause ciliary dyskinesia. [provided by RefSeq, Aug 2015]