

Product datasheet for SC306251

RDM1 (NM_145654) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RDM1 (NM_145654) Human Untagged Clone
Tag:	Tag Free
Symbol:	RDM1
Synonyms:	RAD52B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC306251 representing NM_145654. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGGAGTTGGTACCTTTTGGCGTTCCCATCGAGAGTGACAAAACCTTGCTAGTGTGGGAGCTGAGC
TCCGGACCCACGGCCGAGGCTTTGCATCATTCTCTGTTACAGCATTTTCTCAGTTTGGCCTTCTGTAT
TCAGTCCGGTCTTCCCAAATGCTGCAGTGGCCCATCCTGGTTTCTATGCCGTCATTAAGTTTTATTCT
GCAAGGGTGTCCACAGAGCCCAAAGGCATGCGACCGGAAGCAGCTTTTTCAGAAATCTCCAGTCAAG
GTTTCGTCTTGGCACCAGACATAAGGCAGTTCAACATCAAGCCCTTGCCCTGAACAGTTCCAAATGCCAA
GAACTGGCGAATTACTACTTTGGTTTCAATGGGTGTTCCAAAAGGATCATCAAGTTTCAGGAGCTTTCT
GACCTTGAAGAAAGGGAAAATGAAGATAGCATGGTGCCACTTCCGAAGCAAAGCCTGAAGTTCTTCTGT
GCTTTAGAAAGTGGTGTGCCATCCTGTGATTGCAGGAGTCTGGCATTGGCTTGGTGGAGGAGCCTATG
GATAAGGTGGAGGAAGGACCATTATCATTCTTATGAAAAGGAAGACAGCCAGAAAGCTTGCTATTTCAG
AAGGCTTTGTGATGATGATTCCAGAACTGTTGATTGTTGTTCTAGAAAGTGGTAAAATAGCTGTGGAG
TACAGACCCAGTGAAGACATCGTAGGTGTCAGATGCGAAGAAGAACTACACGGTTTAATTCAGTCCCT
TGCTCTCCCTGGAAGCAGTATGGCCAAGAGGAGGAAGGGTATCTCTCGGATTTACGTTGGAGGAGGAA
GAGTTCAGGCTGCCAGAACTTGACTAG
ACGCGTACGCGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites:	Sgfl-MluI
Plasmid Map:	<input type="checkbox"/>
ACCN:	NM_145654
Insert Size:	855 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_145654.3</u>
RefSeq Size:	1166 bp
RefSeq ORF:	855 bp
Locus ID:	201299
UniProt ID:	<u>Q8NG50</u>
Cytogenetics:	17q12
MW:	32 kDa
Gene Summary:	<p>This gene encodes a protein involved in the cellular response to cisplatin, a drug commonly used in chemotherapy. The protein encoded by this gene contains two motifs: a motif found in RAD52, a protein that functions in DNA double-strand breaks and homologous recombination, and an RNA recognition motif (RRM) that is not found in RAD52. The RAD52 motif region in RAD52 is important for protein function and may be involved in DNA binding or oligomerization. Alternatively spliced transcript variants encoding different isoforms have been reported. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (1) encodes the longest protein (isoform 1, also known as isoform alpha).</p>