

Product datasheet for **SC306293**

NEK11 (NM_145910) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NEK11 (NM_145910) Human Untagged Clone
Tag:	Tag Free
Symbol:	NEK11
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC306293 representing NM_145910. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTGAAATCCAAGAGGCAGCTAAGTGTGTGAGTGGATCAACAGCCATTTCCACTTATCCAAAGACC
TTGATTGCAAGAAGATACGTGCTTCAACAAAACTTGGCAGTGGAAAGTTTGGAACTGTCTATCTGGTT
TCAGACAAGAAAGCCAAACGAGGAGAGGAATTAAGGTAATCTGTTGGAGAACTAAAT
CCAAATGAAACTGTACAGGCCAATTTGGAAGCCCAACTCCTCTCCAAGCTGGACCACCCAGCCATTGTC
AAGTCCATGCAAGTTTGTGGAGCAAGATAATTTCTGCATTATCACGGAGTACTGTGAGGGCCGAGAT
CTGGACGATAAAATTCAGGAATATAACAAGCTGGAAAAATCTTCCAGAAAAATCAAATAATAGAATGG
TTTATCCAGCTGCTGCTGGGAGTTGACTACATGCATGAGAGGAGGATACTTCATCGAGACTTAAAGTCA
AAGAATGTATTTCTGAAAAATAATCTCCTTAAATTTGGAGATTTTGGAGTTTCTCGACTTCTAATGGGA
TCCTGTGACCTGGCCACAACCTTAACTGGAACCTCCCATATATGAGTCCTGAGGCTCTGAAACACCAA
GGCTATGACACAAAGTCGGACATCTGGTCACTGGCATGCATTTTGTATGAGATGTGCTGCATGAATCAT
GCATTCGCTGGCTCCAATTTCTTATCCATTGTTTTAAAAATTTGTTGAAGGTGACACACCTTCTCCCT
GAGAGATATCCAAAAGAACTAAATGCCATCATGGAAAGCATGTTGAACAAGAATCCTTCATTAAGACCA
TCTGCTATCGAAATTTAAAAATCCCTTACCTTGATGAGCAGCTACAGAACCTAATGTGTAGATATTCA
GAAATGACTCTGGAAGACAAAAATTTGGATTGTCAGAAGGAGGCTGCTCATATAATTAATGCCATGCAA
AAAAGGATCCACCTGCAGACTCTGAGGGCACTGTGAGAAGTACAGAAAAATGACGCCAAGAGAAAGGATG
CGGCTGAGGAAGCTCCAGGCGGCTGATGAGAAAGCCAGGAAGCTGAAAAAGATTGTGGAAGAAAAATAT
GAAGAAAAATAGCAAACGAATGCAAGAAATGAGATCTCGGAACTTTGAGCAGCTGAGTGTGATGACTC
CATGAAAAACACATTTAAAAGGAATGGAAGAAAAAGGAGGAGCAACCTGAGGGAAGACTTCTTGTTC
CCCCAGGACGAGGATGAAGAGAGGTGGCAAGGCAGGGAAGAGGAATCTGATGAACCAACTTTAGAGAAC
CTGCCTGAGTCTCAGCCTATTCTTCCATGGACCTCCAGAACTTGAATCAATTGTAGAGGATGCCACA
TCTGACCTTGGATACCATGCAACTCACAGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



[View online >](#)

Restriction Sites:	Sgfl-Mlul
Plasmid Map:	□
ACCN:	NM_145910
Insert Size:	1413 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_145910.3
RefSeq Size:	2107 bp
RefSeq ORF:	1413 bp
Locus ID:	79858
UniProt ID:	Q8NG66
Cytogenetics:	3q22.1
Protein Families:	Druggable Genome, Protein Kinase
MW:	54 kDa
Gene Summary:	<p>This gene encodes a member of the never in mitosis gene A family of kinases. The encoded protein localizes to the nucleoli, and may function with NEK2A in the S-phase checkpoint. The encoded protein appears to play roles in DNA replication and response to genotoxic stress. Alternatively spliced transcript variants have been described.[provided by RefSeq, Mar 2009]</p> <p>Transcript Variant: This variant (2) includes an alternate exon, compared to variant 1, that causes a frameshift. The resulting protein (isoform 2), also known as NEK11S, is shorter and has a distinct C-terminus compared to isoform 1.</p>