

Product datasheet for **SC306311**

WFDC9 (NM_147198) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	WFDC9 (NM_147198) Human Untagged Clone
Tag:	Tag Free
Symbol:	WFDC9
Synonyms:	dj688G8.2; WAP9
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_147198, the custom clone sequence may differ by one or more nucleotides

ATGAAGCCCTGGATTCTTCTACTCGTCATGTTTCATCTCTGGAGTTGTGATGCTTCTGCCTGTGCTGGGAA
 GCTTCTGGAACAAAGATCCCTTTCTAGATATGATAAGAGAACTGAGCAGTGCTGGGTACAGCCTCCATA
 TAAGTACTGTGAGAAAAGGTGTACTAAAATAATGACTTGTGTACGTCCAAATCATACATGCTGCTGGACC
 TACTGTGGAAACATCTGCTTAGACAACGAAGAGCCCTTAAATCAATGCTAAACCCCTAG

Restriction Sites:	Please inquire
ACCN:	NM_147198
OTI Disclaimer:	Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)


[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_147198.2 , NP_671731.1
RefSeq Size:	659 bp
RefSeq ORF:	270 bp
Locus ID:	259240
UniProt ID:	Q8NEX5
Cytogenetics:	20q13.12
Protein Families:	Secreted Protein, Transmembrane
Gene Summary:	The WAP-type four-disulfide core (WFDC) domain, or WAP signature motif, contains eight cysteines forming four disulfide bonds at the core of the protein, and functions as a protease inhibitor in many members of the WFDC domain family. This gene encodes a protein which contains a WFDC domain, and is thus a member of the WFDC domain family. This gene and several other gene family members are clustered at 20q13.12. [provided by RefSeq, Jul 2008]