

Product datasheet for **SC306373**

ZNF441 (NM_152355) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF441 (NM_152355) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF441
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF:	>NCBI ORF sequence for NM_152355, the custom clone sequence may differ by one or more nucleotides <pre> ATGGACTCAGTGGCATTGAGGATGTGGCTATAAACTTCACCTGTGAGGAGTGGGCTTTG CTGGGTCCATCACAGAAGAGTCTCTACAGAGATGTGATGCAGGAAACCATCAGAAACCTG GACTGTATAGGAATGATATGGCAAAATCATGATATAGAAGAAGATCAGTACAAAGATCTC AGAAGAAATCTAAGATGTCATATGGTAGAGAGAGCCTGTGAAATTAAGATAATAGTCAA TGTGGAGGACCCTTTACCCAGACTCAAGACAGTATTGTGAACGAGAAAAATACCTGGAGTA GATCCATGGGAAAGCAGTGAGTGTACAGACGTCCTCATGGGTCGTTTCATCTCTTAATTGC TACGTTAGAGTTGACAGTGAACACAAACCATGTGAGTATCAAGAATATGGAGAGAAGCCA TATACACATACACAATGTGGGACAGCTTTCAGTTATCAGCCCTGCTTCAAATACATGAA AGACCTCAGCATGGAAAGAACTCTATGATTGTAAGGAATGTGCAAGCTTCAGTTCTCTT GAAAACCTTCAAAGACACATGGCAGCACACCATGGAGATGGACCTCGTATATGTAAGTTG TGTGGAAACGCCTTTATTTGGCCTAGTTTATTTTCATATGCTTAGAAGAACTCACACTGAA GAGAAACCATATGAATATGAGCAGTGTCTACAGCGTTTCCTGCTTATAGTTCCACTCTA AGACATGAAAGAACACACAGTGGAGAGAAACCTATCAATGTAACAATGTGGGAAAGCC TTCAGTTGTTCTGTTTACACTCAACTATATGAAAGGACTCACACTGGAGAACAACTCTAT GAATGTAAGCAATGTGGGAAAGCATTATCATCTTGAAGCTTTCAAAGACACATGATA GTGCACACTGGAGATGGGCCTCATAAATGTAAGATATGTGGAAAAGGCTTCTTTCTCCC AGTTTCAGTTTGAAGACATAAAAGAACTCACACTGGAGAGAAACCGTATGAATGTAAGTAT TGTGGGAAAGCATTCTCTGATTGCACAGGTTTTCGAAGACACATGATAACGCACACTGGA GATGGACCTCATAAATGCAAGGTATGTGGGAAAGCCTTTGATTCTCCAGTTTATGTCGA AGGCATGAAACAACTCATACTGGGGAGAAACCTATAAATGTGAATGTGGGAAAGCCTTT AGTGATTTCTATTACTTTGAAATCATGAACTACTCACACTGGAGAGAAAGCCATATAAA TGTAACAATGTGGAAAAGCCTTCATTTGTTGCACTTACCTTCAAATACATGAAAGAATT CACACTGGGGAGAGACCCTATAAATGTAACAATGCGGAAAAGCCTTCAGGTCTTCCAAT TACATTCGAGTACATGAAAAGACTCACACTGGAGAAAAGCCTATGAATGTAAGCAGTGT GGAAAAGCACTTCTCATCTGAAAAGCTTTCAGAGACACATGATAATGCACACTGGAGAT GGACCTCATAAATGTAAGATATGTGGGAAAAGCTTTCAGTTCTCCAGTTTCATTTCGAAGA CATGAAAGAATTCACACTGGGGAGAGACCCTATAAGTGTAACTATGTGGGAAAGCCTTC AGGTCTTCCAGTTACATTCAACTACATGAAAGGACTCACACTGGAGAGAAACCTATGGT TGTCAGCAATGTGGGAAAGCATTATCTGATCTCTCAAGCTTTCGAAGACACATGATAACA CATACTGGAAATGGACCTCATAAATGTAAGATATGTGGGAAAAGGCTTTCAGTTATCCAGT TCAGTGCAAAGACATGAAAGAACTCACACTGGAGAGAAACCTATGAATGCAAGGAATGT GGTAAAGCCTTCAGTCATTCAAGTTACTTACGAATACACGAAAGAGTTTCACTGGAGAG AAGCCGTATAAATGTAAGGAATGTGGGAAACCATTCATTGTCCAGTGCCTTTCATAAA CATGAAAGGACCCACAGTATGGAGAAACCTATAAGTGTAAAGAATGTGGGGAAGCATT CATTGTATCAGTTCCTTTCATAAACATGAAATGACTCACTAG </pre>
Restriction Sites:	Please inquire
ACCN:	NM_152355
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_152355.1, NP_689568.1</u>
RefSeq Size:	2829 bp
RefSeq ORF:	1881 bp
Locus ID:	126068
UniProt ID:	<u>Q8N8Z8</u>
Cytogenetics:	19p13.2
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]