

Product datasheet for SC306399

RAD9 homolog B (RAD9B) (NM_152442) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	RAD9 homolog B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PSI00001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC306399 representing NM_152442. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAAATTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCAGCCATGCTGAAGTGCCTGATGAGCGGCAGTCAGGTGAAAGTATTTGGGAAAGCAGTTCAAGCT
CTATCACGAATTAGTGACGAGTTCTGGCTAGACCCATCTAAAAAGGTCTTGCTCTAAGATGTGTGAAT
TCTTCTCGGTGAGCATATGGATGTGCTGTTCTCTCTGTGTTTTTCAGCATTATCAATGGTCAGCT
TTAGTGAAATGAGTGAAATGAACCTGACACAACACTGCATTTAAATGCAAATTGGGAATGAAGTCA
ATTTTGCCCATCTTTAGATGTCTGAATCCCTTGAAAGAAATATAGAGAAGTGCAGAATATTCACCAGA
TCTGATAAATGCAAAGTAGTTATTCAATCTTCTACAGACATGGTATTAAGAAGCTCATAATATATGT
TTTCAAGAAAGTCAGCCTTTGCAAGTTATTTTGACAAGAATGTTGTACTAATACGCTAATGATTCAA
CCAAGATTGCTTGCTGATGCCATTGTTCTTTTACATCAAGTCAAGAGGAAGTTACTCTTGCTGTTACT
CCACTGAATTTTGCCTCAAGAGTTCTAATGAGGAATCAATGGATTTGAGCAATGCTGTACACAGTGAG
ATGTTTGTGGCTCAGATGAGTTTGACTTCTTCAAATTGGAATGGACACTGAGATAACATTTTGTTC
AAAGAATTGAAGGAATACTGACATTTTCAGAAGCTACACATGCTCCTATATCCATTTATTTGATTTC
CCTGGGAAACCTCTGGCTTTGAGTATTGATGATATGTTAGTGGAAGCTAACTTTATTTTGGCCACATTA
GCTGATGAACAAAGTAGAGCATCTTACCACAGTCAGTGTGCTTTTACAGAAACGAAAAAGGTCAGAT
CTGATTGAAAAAAGGCTGGCAAAATGTAAGTGGCCAGGCCCTGGAATGTATTTCAAAAAAGCAGCA
CCAAGAAGGCTTTATCTAAGGAGACTCTCACAACATATCTGCATTGGAAAAGTGTGGCAGCCCTGCA
ATGAAAAGAGTGGATGGAGATGTCAGTGAAGTATCAGAAAGCAGTGTGAGCAACACAGAGGAAGTGCCA
GGGTCTCTGTGCTCAGAAAGTTTTCTGTCATGTTCTTTGGAGCAGTTTCTTCTGACCAGCAAGAACAC
TTCAACCACCCTTTGACAGTCTGGCAAGAGCAAGTGACAGTGAAGAGGACATGAATAATGTGTGCTGC
AGGAAAGAATTTAATGGAAGTATGCCAAATTTCTGTATTATCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



ACCN:	NM_152442
Insert Size:	1290 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_152442.3
RefSeq Size:	1898 bp
RefSeq ORF:	1290 bp
Locus ID:	144715
UniProt ID:	Q6WBX8
Cytogenetics:	12q24.11
MW:	48.1 kDa