

Product datasheet for **SC306404**

C17orf77 (NM_152460) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C17orf77 (NM_152460) Human Untagged Clone
Tag:	Tag Free
Symbol:	C17orf77
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC306404 representing NM_152460. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
 ATGGATGAATTGGCGCTGTCATTTTCCCTGACATGTCTCCTTCTGAGAACAGAGCGTCTGTCTCCA
 AGCCAGCCGCTGTCAATCCAGTGCCTTAAAGCTCCAGCCACACTCACCTGGGAGGATGAGAAGCAGCAG
 AGATGGGACAGCCACATGGTCCTGTCTCCTCTCCTTGGTGATCAGAGTGTCTGGTGCCCTTC
 AGGGACTTGAATCCAAGCTCGGAAGTCAACACCGCAAACCTACTCGAATCTCCTTCAAGTCTCTTGCTC
 ACTTCTGTGTACATCTGTAGCTACTTCTTTTTACATATTAGGAGAAAAGAGGTGTCTTCTTTGAAA
 AGACTAAGGTACAGTGTATGTTGCAAGGTGTGCCAAACTTTTGTGCATGTGAAAAGAAAATGTGTCC
 GGGACTGGTCAGGTTTGCACGGGAGTGCACGTAGGAGCCAAGGAGCAGGAGGAGCCAGGAGGAACCCAA
 GCATTGAGGAGCTGCGGGATCTACTGCCTGGAGGAGAGGACGGACAAGGCCAGCCATGAGGAGTGCAGG
 GAGAGGAGCACTTTGGGAAGACCACAGTGCACAGGACTCACACCCAGCCTTGCTGGAGAAAGCCCATGC
 CCCAGGCTTCTCCAGGCAGCCCACTGTAAGACATCTCATTGCCTTCTCCTGCCCTGGGCTGTCTGAC
 CCACTGCCCTTCCCCAGGAACACTTCCACTGGGATCA**TGA**
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_152460
Insert Size:	732 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_152460.2
RefSeq Size:	2525 bp
RefSeq ORF:	732 bp
Locus ID:	146723
UniProt ID:	Q96MU5
Cytogenetics:	17q25.1
Protein Families:	Secreted Protein
MW:	26.3 kDa