

Product datasheet for **SC306406**

SLC35G3 (NM_152462) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC35G3 (NM_152462) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLC35G3
Synonyms:	AMAC1; TMEM21A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC306406 representing NM_152462. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGCATCGCC**
 ATGGCTGGCAGTCACCCCTATTTCAACCAGCCTGACTCCACACACCCATCGCCGCCCTCCGCTCCACCC
 AGCCTCCGCTGGTACCAGCGCTGCCAGCCCTCTGATGCCACCACTGGCCTGCTGGTGGCCCTGCTGGGT
 GGGGGCCTGCCTGCTGGCTTCGTGGGCCCCCTTCTCGTATGGCTTACCAGGCTTCCAACCTGCCCTCG
 CTGGAGCTGCTCATCTGGCGATGCCTCTCCACCTCCCTATTGCCCTGCTACTTAACTGCGTGGCGAC
 CCCCTTCTGGGAACCTCTGACATCCGAAGCCGGGCCTTCTCTGTGCCCTGCTCAACATCCTCAGCATT
 GGATGTGCCTACAGTGCCTTCAGGTGGTGGCCGCTGGCAACGCTGCCACTGTTCCGAAAGTTCTTCC
 ACCGTCTGCTCCGCCGTCTCACTCTCTGCCTTGAGAGCCAGGGTCTCAGTGGCTACGACTGGTGTGGA
 CTGTTGGGCTGCATCCTAGGACTAATCATCATTGTGGGACCTGGACTCTGGACACTACAGGAGGGGACC
 ACGGGTGTCTACACCGCCCTGGGCTATGTGGAGGCTTCTCTGGGAGGCCTGGCGCTGTCCCTGAGGCTT
 CTGGTCTATCGTTCTCTGCACTTCCCCCTGCCTCCCAACAGTGGCCTTCTATCTGGCTTGGTGGG
 CTGCTGGGCTCTGTGCCAGGCCTCTTGTGCTGCAGGCCCGCTGTTGCCAGTGACCTCCTGAGTTGG
 AGTTGTGTGGGGCAGTGGGGATCCTCGCTTGGTCTCCTTCACATGTGTGGGCTATGCGGTACCAAG
 GCCACCTGCCCTGGTGTGCGCTGTCTACATTCCGAGGTGGTGTGGCCCTTATACTGCAGTATTAT
 ATGCTCCATGAGACTGTGCACCTTCTGACATCGTGGCGCAGGGGTTGTGCTGGGCAGCATTGCCATC
 ATTACAGCCCAGAACCTCAGCTGTGAGAGGACAGGGAGGGTGGAGGAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_152462
Insert Size:	1017 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_152462.2</u>
RefSeq Size:	1874 bp
RefSeq ORF:	1017 bp
Locus ID:	146861
UniProt ID:	<u>Q8N808</u>
Cytogenetics:	17q12
Protein Families:	Transmembrane
MW:	35.4 kDa