

Product datasheet for **SC306414**

SAMD11 (NM_152486) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SAMD11 (NM_152486) Human Untagged Clone
Tag:	Tag Free
Symbol:	SAMD11
Synonyms:	MRS
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC306414 representing NM_152486.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGATCGCC
ATGTCCAAGGGGATCCTGCAGGTGCATCCTCCGATCTGCGACTGCCCGGGCTGCCGAATATCCTCCCCG
GTGAACCGGGGGCGGCTGGCAGACAAGAGGACAGTCGCCCTGCCCGCCCGAACCTGAAGAAGGAG
CGAACTCCAGCTTCTCTGCCAGCGATGGTGACAGCGACGGAGTGGCCCCACCTGTGGCGGCGGCCA
GGCTTGAAGCAGGAGGATGGTCCGCACATCCGTATCATGAAGAGAAGAGTCCACACCCACTGGGACGTG
AACATCTCTTTCCGAGAGGCGTCTGCAGCCAGGACGGCAACCTTCCACCCCTCATATCCAGCGTCCAC
CGCAGCCGCCACCTCGTTATGCCCGAGCATCAGAGCCGCTGTGAATTCAGAGAGGCAGCCTGGAGATT
GGCCTGCGACCCGCCGTGACCTGTTGGCAAGAGGCTGGGCCGCTCCCCCGTATCAGCAGCGACTGC
TTTTCAGAGAAGAGGGCACGAAGCAATCGCTCAAGAGGCGCTGCTGCTGCCCGGGAGCTGGGGCCC
AGCATGGCCCCGGAGGACCATTACCGCCGGCTTGTGTAGCACTGAGCGAGGCCAGCACCTTTGAGGAC
CCTCAGCGCTCTACCACCTGGGCTTCCAGCCACGGTGAGGACCCACCCTGGCATGATCCCCCTCAT
CACCTCCCAGCCAGCATCTCTGAGGGTCCGGCAGGAGGTGGCGGCTGCAGCTCTGAGGGGCCCACT
GGCCTGGAAGCCCACTGCCCTCTCCACGGCAGGTGAGCGTCGGAAGCAGGGCTGGCTCAGCACCCGG
GAGGGCGCCGCCAGCTGCCGCCCGTCTTCTCGGAGAGGGAGCTGCCCTCAGCCGCCCCCTTGTCTG
TCGCCGAGAATGCCCTCACGTGCCCTGGGCCCCATCTCAGGCCCCCTTCTGGGGTGCCCTCG
GCTCTGTGCCAGACCCAGGCTACGGCTTCTGCCCGCCGCGAGGCGGAGATGTTGCCCGGCAGCAG
GAGCTCCTGCGGAAGCAGAACCTGGCCCGGCTGGAGCTGCCCGCCGACCTCCTGCGGCAGAAGGAGCTG
GAGAGCGCGGCCACAGCTGCTGGCGCCGAGACCGCCTGCGCCCCAACGACGGCGCCGAGGAGCTG
CAGCGGCGCGGGCCCTGCTGGTGTGAACACGGCGCGGCCACTGCTGGCCCTGCCCCCAGGGG
CCCCGGGCTCCGAGCCCCCACCCTGTCGGGACTCTGCCCGCGAGCCCCCGGAAGGGGGTCCC
GGCCCTGCCTCAGCGCGGCCAGCGAGTCCAAGGAGATGACGGGGCTAGGCTCTGGGCACAAGATGGC
TCGGAAGACGAGCCCCCAAAGACTCGGACGGAGAGGACCCGAGACGGCAGCTGTTGGGTGACGGGG
CCCACTCCGGGCAAGCTCCAGCTGGAGGGGCCGCGCCGAGGGGAAGGGGCTTTTCCAGGGTCCACA
CTGCCCTGGGCTTCCCTATGCCGTGAGCCCTACTTCCACACAGGCGCGGTAGGGGACTCTCCATG
GATGGGGAGGAGGCCCGACCCCTGAGGACGTACCAAGTGACCGTGGATGACGTCTGCAGTTCGTG
GGGGCCCTGTCTGGCTGTGAGAGTAACTCGGGTCTTACGGGAGCAGGGGATCGACGGGAGACCTG
CCACTGCTGACGAGGAGCACCTGCTGACCAACATGGGGCTGAAGCTGGGGCCCGCCTCAAGATCCGG
GCCAGGTGGCCAGGCGCTGGGCCGAGTTTCTACGTGCCAGCTTCCCGTGGCTCTGCCACTGCAG
CCACCAACCCTGCGGGCCCCGGAGCGAGAACTCGGCACAGGAGAGCAGCCCTGTCCCCCAGCAGGCC
ACGTCCCCCTATGGAGGGGGCCACGCCCTTGCCGGTCAAACCTTACCCAAAGCAGGAGAATGGGACCTTG
GCTCTACTTCCAGGGGCCCGACCTTCCAGCCTCTGTGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_152486

Insert Size: 2046 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_152486.2</u>
RefSeq Size:	2575 bp
RefSeq ORF:	2046 bp
Locus ID:	148398
UniProt ID:	<u>Q96NU1</u>
Cytogenetics:	1p36.33
MW:	72.7 kDa
Gene Summary:	May play a role in photoreceptor development.[UniProtKB/Swiss-Prot Function]