

Product datasheet for **SC306416**

UBE2U (NM_152489) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: UBE2U (NM_152489) Human Untagged Clone
Tag: Tag Free
Symbol: UBE2U
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >SC306416 representing NM_152489.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCACGGCAGAGCTTACCTCTTGCTGCACAGAGACTTCTGTGATCTCAAGGAGAACAATTATAAGGGT
ATCACTGCTAAGCCTGTAAGTGAAGATATGATGGAATGGGAAGTTGAAATTGAAGGTCTACAGAATTC
GTTTGGCAGGGTTTAGTCTTCCAACGACAATACATTTTACATCGGAGTACAACATGCTCCTCCAGTT
GTGAAATTTATAACAATCCGTTTCATCCAATGTAGACCCACACACTGGTCAGCCCTGTATAGACTTT
TTGGACAACCCCTGAGAAGTGAATACAACTATACATTGAGCAGCATCTTACTTGCCTACAGGTTATG
CTTTCTAATCCAGTGCTAGAGAATCCAGTGAATTTGGAAGCAGCCAGAATACTGGTTAAAGATGAATCT
CTGTACAGAACAATTCTAAGACTTTTCAACAGGCCATTACAAATGAAAGATGACAGCCAGGAGTTACCT
AAAGACCCACGTAATGTATCAGACCAATTAACAACCTCATTTAGTGATTACTACCAGACATGGTCC
AGAATAGCTACATCAAAGCCACAGAATACTACAGAACTCCATTGCTCAAAGTTCCAAATTTTCATTGGA
CAGTATTACAAATGGAAGAAAATGGATCTACAGCATCAGAAAGAATGGAATTTAAATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map: □

ACCN: NM_152489

Insert Size: 681 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_152489.1
RefSeq Size:	1163 bp
RefSeq ORF:	681 bp
Locus ID:	148581
UniProt ID:	Q5VWX9
Cytogenetics:	1p31.3
Protein Pathways:	Ubiquitin mediated proteolysis
MW:	26.7 kDa
Gene Summary:	Catalyzes the covalent attachment of ubiquitin to other proteins.[UniProtKB/Swiss-Prot Function]