

Product datasheet for **SC306454**

PGBD4 (NM_152595) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PGBD4 (NM_152595) Human Untagged Clone
Tag:	Tag Free
Symbol:	PGBD4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC306454 representing NM_152595.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGTCAAATCCTAGAAAACGTAGCATTCTATGCGTGATAGTAATACCGGTCTCGAACAGTTGTTGGCT
GAAGATTCATTTGATGAATCTGATTTTTCGAAATAGATGATTCTGATAATTTTTCGGATAGTCTTTA
GAAGCCGATAAGATCAGGCCTCTGTCCATTTAGAATCTGATGGAAAGAGCTCTACATCAAGTGACTCA
GGGCGCTCCATGAAATGGTCAGCTCGTGCTATGATTCCACGTCAAAGGTATGACTTTACCGGCACACCT
GGCAGAAAAGTCGATGTCAGTGATATCACTGACCCATTGCAGTATTTGAACTGTTCTTTACTGAGGAA
TTAGTTTCAAAAATTACTAGAGAAACAAATGCCAAGCTGCCTTGTGGCTTCAAAGCCACCGGGTCCG
AAAGGATTTTCGGAATGGATAAATGGAAGACACTGACAATGACGAGCTCAAAGCTTTTTTGCAGTA
ATGTTACTGCAAGGTATTGTGCAGAAACCTGAGCTGGAGATGTTTGGTCAACAAGGCCTCTTTGGAT
ACACCTTATCTCAGGCAAATTATGACTGGTGAAGATTTTACTTTTGTTCGGTGCCTGCATTTTGTG
AACAACTCTTCTATATCTGCTGCTCAATCAAAGGCCAGATTTTCATTGCAGAAGTCAAACCTGTGTTC
GACTTTCTGTAAATAAATTTTCCACTGTATATACTCCAAACAGAAACATTGCAGTTGATGAATCACTG
ATGCTGTTCAAGGGGCCATTAGCTATGAAGCAGTACCTCCCGACAAAACGAGTACGATTTGGTCTGAAG
CTATATGTAATTTGTGAAAGTCAGTCTGGTTATGTGTGGAATGCGCTTGTTCACACAGGGCCTGGCATG
AATTTGAAAGATTCAGCGGATGGCCTGAAATCATCACGCATTGTTCTTACCTTGGTCAATGACCTTCTT
GGCCAAGGGTATTGTGTCTTCTCGATAACTTTAATATATCTCCCATGCTTTTCAGAGAATTACATCAA
AATAGGACTGATGCAGTTGGGACAGCTCGTTTGAACAGAAAACAGATTCCAAATGATCTGAAAAAAGG
ATTGCAAAGGGGACGACTGTAGCCAGATTCTGTGGTGAAGTTATGGCACTGAAATGGTGTGACGGCAAG
GAGGTGACAATGTTGTCAACATTCACAATGATACTGTGATTGAAGTAACAATAGAAATGGAAGAAA
ACTAAAAGGCCACGTGTCAATTGTGGATTATAACGAGAATATGGGAGCAGTGGACTCGGCTGATCAATG
CTTACTTCTTATCCATCTGAGCGCAAAAGACACAAGGTTTGGTATAAGAAATCTTTACCATCTTCTA
CACATTACAGTGCTGAAGTCTTACATCTGTTCAAGAAGGATAATCCTGAGCACACGATGAGCCATATA
AACTTCAGACTGGCATTGATTGAAAGAATGCTGGAAAAGCATACAAGCCAGGGCAGCAACATCTTCGA
GGTCGTCTTGTCCGATGATGTACACCTCTTCGTCTGTCTGGAAGACATTTCCCAAGAGCATACCA
GCAACGTCGGGAAACAGAATCCAATGGTCGCTGCAAAATTTGCTGCTCCCAATACGACAAGGATGGC
AAGAAGATCCGGAAGAAACGCGCTATTTTGTGCCGAATGTGATGTTCCGCTTTGTGTTGTTCCGTGC
TTTGAAATTTACCACACGAAAAAAATTATA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_152595

Insert Size: 1758 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_152595.4](#)

RefSeq Size: 2318 bp

RefSeq ORF: 1758 bp

Locus ID: 161779

UniProt ID: [Q96DM1](#)

Cytogenetics: 15q14

MW: 67 kDa

Gene Summary: The piggyBac family of proteins, found in diverse animals, are transposases related to the transposase of the canonical piggyBac transposon from the moth, *Trichoplusia ni*. This family also includes genes in several genomes, including human, that appear to have been derived from the piggyBac transposons. This gene belongs to the subfamily of piggyBac transposable element derived (PGBD) genes. The PGBD proteins appear to be novel, with no obvious relationship to other transposases, or other known protein families. [provided by RefSeq, Jul 2008]