

Product datasheet for **SC306470**

CCER1 (NM_152638) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CCER1 (NM_152638) Human Untagged Clone
Tag:	Tag Free
Symbol:	CCER1
Synonyms:	C12orf12
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC306470 representing NM_152638. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGCATCGCC**
 ATGACTCAGACCTCGACACAAGGGAAGACCTCTGAACCTGGGCGGCGGTGGCGCGCGGCTGTGGC
 TGTGGCTGGGCACACTCGGCCTCCTTGAGCTCCTGGTCGTCCTGCCATCGAAGGCGCCGGGTGCTCCA
 GCGTACAATAGACCGCACCGATATAGCCCCAAGACCGAGTATGGGCCCCAAGGAAGCAGCCGAAGCAA
 CAGCACGGCCCGGGCTTTTGGTTCCAACCCCGTGTGCTCTAACTGGGGGTGCTGGGGAGGGCCCTGG
 CGCCACCCCTCCAGGATCTGGAAATCCCCTGCCCGGTGCAAGTGTTCGGGTGTATGGCCTGCAC
 CCTCTGCTTTTGTGCTGCTCCTGCTGGAGCGGGTCTGGAACCTGGCTGGTGAAGCCCCAGGC
 AGGAAGAAGCGCTGGGGCCGAGAGGCCCGGCCCTGCCACCACCTCGCCACTCCTACCCGCGGAGC
 CCGCCAGCGGATGTGAGCACGCTGCCCGCGCGGTCAAGCTGTATGAGTGGAGAGAGCCTGGCATGCGA
 GCGCCGCCCAACACCACCAATTCATCATGAACCAGATCTACGAGGACATGAGGCAGCAGGAGAAGGTG
 GAGCGTCAGCAGGAGGCGCTGCGGGCGCAGAAGGCCACGGTGAGCGGCGAGGCTCCCCAGCCAGATCC
 TCCGAAACGACGCGCCCTGGCGGCAGCAAGAAACCTGGGGACTGCAGGAACTCTGTATGGCTTT
 GTGCAGAATCCCTCTAGCATTAGTCCCAACCCAGAGGAAAACCAAGTCTTTGCCCGCTGCTGGTG
 GAAGAAGAGGAGGAGAAGAAAAATGATGATGAGGAGGAGTATGACCAGGAGGTGTGTATGCAAAGGAG
 GCGAGCGAGGAGGAAGAAGAGGTCAAGATGAGGAGGAAGAGGTCAAGATGAGGAGGAAGAAGAGGTC
 GAAGAGGCTGAATATGTGGAGGAGGAGAGGAGGAGCTGGAAGAGGAGGAGCTGGAAGAGGAAGAGGAG
 GTCTTGAGGAGAACGAGCAGAGAGGGGAAGAATTTCACTTGCCTCTGGAATGCCTTTATCAATCTTC
 GTAGAGGCTGAAGAAAAGAGAGAGAATTTATAAGCTGCACTTTTTTAAACCCAGAGCAGATAATTCCC
 AAAGTGCCACAGGAATCCCTGTTTCATGGCACAGGACTTTAACTGT**TAG**
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC

Restriction Sites: SgfI-MluI


[View online »](#)

ACCN:	NM_152638
Insert Size:	1221 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_152638.3</u>
RefSeq Size:	2977 bp
RefSeq ORF:	1221 bp
Locus ID:	196477
UniProt ID:	<u>Q8TC90</u>
Cytogenetics:	12q21.33
MW:	46.5 kDa