

Product datasheet for **SC306496**

CNPY4 (NM_152755) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CNPY4 (NM_152755) Human Untagged Clone
Tag:	Tag Free
Symbol:	CNPY4
Synonyms:	PRAT4B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC306496 representing NM_152755. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGACCTGTGCGGTTGGGAATATTGCTTTTCTTTTTTTGGCCGTGCACGAGGCTTGGGCTGGGATG
TTGAAGGAGGAGGACGATGACACAGAACGCTTGCCAGCAAAATGCGAAGTGTGTAAGCTGCTGAGCACA
GAGCTACAGGCGGAACTGAGTCGCACCGGTCGATCTCGAGAGGTGCTGGAGCTGGGGCAGGTGCTGGAT
ACAGGCAAGAGGAAGAGACACGTGCCTTACAGCGTTTCAGAGACAAGGCTGGAAGAGGCCTTAGAGAAT
TTATGTGAGCGGATCCTGGACTATAGTGTTCACGCTGAGCGCAAGGGCTCACTGAGATATGCCAAGGT
CAGAGTCAGACCATGGCAACACTGAAAGGCCTAGTGCAGAAGGGGTGAAGGTGGATCTGGGGATCCCT
CTGGAGCTTTGGGATGAGCCAGCGTGGAGGTCACATACCTCAAGAAGCAGTGTGAGACCATGTTGGAG
GAGTTTGAAGACATTGTGGGAGACTGGTACTTCCACCATCAGGAGCAGCCCTACAAAATTTTCTCTGT
GAAGGTGATGTGCTCCCAGCTGCTGAACTGCATGTCTACAGGAACTTGGACTGGAAAGGAGATCACA
GATGGGGAAGAGAAAACAGAAGGGGAGGAAGAGCAGGAGGAGGAGGAAGAGGAGGAAGAGGAAGGG
GGAGACAAGATGACCAAGACAGGAAGCCACCCAACTTGACCGAGAAGATCTTTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites:	Sgfl-Mlul
Plasmid Map:	<input type="checkbox"/>
ACCN:	NM_152755
Insert Size:	747 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_152755.1
RefSeq Size:	1547 bp
RefSeq ORF:	747 bp
Locus ID:	245812
UniProt ID:	Q8N129
Cytogenetics:	7q22.1
MW:	28.3 kDa
Gene Summary:	Plays a role in the regulation of the cell surface expression of TLR4.[UniProtKB/Swiss-Prot Function]