

## Product datasheet for **SC306584**

### DEFB119 (NM\_153323) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	DEFB119 (NM_153323) Human Untagged Clone
Tag:	Tag Free
Symbol:	DEFB119
Synonyms:	DEFB-19; DEFB-20; DEFB20; DEFB120; ESC42-RELA; ESC42-RELB
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_153323, the custom clone sequence may differ by one or more nucleotides ATGAAACTTCTTTACCTGTTTCTTGCCATCCTTCTGGCCATAGAAGAACCAGTGATATCA GTAGAGTGTGGATGGATGGACTGCCGGTGTGTGCAAAGATGGTGAAGACAGCATC ATACGCTGCCGAAATCGTAAACGGTGCTGTTCCTAGTCGTTATTTAACAATCCAACCA GTAACAATTCATGGAATCCTTGGCTGGACCACTCCTCAGATGTCCACAACAGCTCCAAAA ATGAAGACAAATATAACTAATAGATAG
Restriction Sites:	Please inquire
ACCN:	NM_153323
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>



[View online »](#)

RefSeq: [NM\\_153323.3](#), [NP\\_697018.1](#)

RefSeq Size: 467 bp

RefSeq ORF: 267 bp

Locus ID: 245932

UniProt ID: [Q8N690](#)

Cytogenetics: 20q11.21

Protein Families: Secreted Protein

**Gene Summary:** This gene encodes a member of the beta subfamily of defensins. Beta-defensins are antimicrobial peptides that protect tissues and organs from infection by a variety of microorganisms. This gene is found in a cluster with other beta-defensin genes on the long arm of chromosome 20. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2014]  
Transcript Variant: This variant (3) encodes the longest isoform (c). Isoform c is thought to be an antimicrobial peptide.