

Product datasheet for **SC306647**

C9orf96 (STKLD1) (NM_153710) Human Untagged Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | C9orf96 (STKLD1) (NM_153710) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | STKLD1 |
| Synonyms: | C9orf96; SgK071; Sk521 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |



[View online »](#)

Fully Sequenced ORF: >SC306647 representing NM_153710.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGCTTTGGGCCAGGGTCCAATCGCAGGGCCACGCAGGGGAGCGAGGCCAGGGTCCCCGGAGAG
CCCATGGAGAAGTACCAGGTTTTGTACCAGCTGAATCCTGGGGCCTTGGGGTGAACCTGGTGGGAG
GAAATGGAAACCAAAGTCAAGCATGTGATAAAGCAGGTGGAATGCATGGATGACCATTACGCCAGTCAG
GCCCTGGAGGAGCTGATGCCACTGCTGAAGCTGCGGCACGCCACATCTCTGTGTACCAGGAGCTGTTC
ATCACGTGGAATGGGAGATCTTCTCTGTACCTCTGCCTGGTGTGGAGTTCAATGAGCTCAGCTTC
CAGGAGGTCATTGAGGATAAGAGGAAGGCAAAGAAAATCATTGACTCTGAGTGGATGCAGAATGTGCTG
GGCCAGGTGCTGGACGCGCTGGAATACCTGCACCATTTGGACATCATCCACAGGAATCTCAAACCTCC
AACATCATCTCATCAGCAGTGACCACTGCAAACCTGCAGGACCTGAGTTCCAATGTGCTAATGACAGAC
AAAGCAAATGGAATATTCGTGCGGAGGAAGACCCCTTTCGTAACTCCTGGATGGCCCTGAAGCCCTC
AATTCTCCTTCAGCCAGAAATCAGACATCTGGTCCCTGGGCTGCATCATTCTGGACATGACCAGCTGC
TCCTTCATGGATGGCACAGAAGCCATGCATCTGCGGAAGTCCCTCGCCAGAGCCAGGCCAGCCTGAAG
GCCGCTCTGAAGACAATGGAGGAGAAGCAGATCCCGGATGTGGAACCTTCAGGAATCTTCTGCCCTTG
ATGCTCCAGATCGACCCCTCGGATCGAATAACGATAAAGGACGTGGTGCACATCACCTTCTTGAGAGGC
TCCTTCAAGTCTCGTGCCTCTCTGACCCTGCACCGGCAGATGGTGCCTGCCTCCATCACCGACATG
CTGTTAGAAGGCAACGTGGCCAGCATTTAGAGGTCATGCAGAAATCTCTGGCTGGCCCGAAGTCCAG
CTCAGGGCCATGAAGAGGCTTCTGAAAATGCCTGCAGATCAGTAGGTCTGCCGTGGCCCCGGAGCTG
GTGGAGGTGGTGGTACGACCATGGAGCTACATGACAGGTCCTCGATGTCCAGCTGTGTGCTGTCTCC
CTGTCTGCACCTCCTGGGCCAAGCGTGGTGCACCACCGGAAGCCAAGGCTCCCTGCAACCAAGCC
ATCACCTCCACCCTGCTGAGTGCTCTTCAGAGCCACCCGAGGAGGAGCCACTTCTTGCATGGTCTAC
AGCCTGTAGCCATCACACAACCCAGGAGTCAGAGTCACTGTGAGAGGAGCTGCAGAATGCTGGGCTG
CTGGAGCACATCTGGAGCACCTCAACAGCTCCCTCGAAAGCAGGGACGTCTGCGCCAGCGGCCTGGGC
CTGCTCTGGGCCCTCCTGCTGGACGGTATCATTGTGAACAAGGCCCTTGGAGAAGTCCCGGACCTC
ATCAGCCAGGTGTTGGCCACCTACCCTGCGGATGGGAAATGGCAGAAGCCAGCTGCGGAGTCTTCTGG
CTGCTGTCCCTGTGGGCTGCATCAAGGAGCAGCAGTTTGAACAAGTGGTGGCGTGTCTCTGCAAAGC
ATCCGGCTGTGCCAGGACAGAGCCCTGCTGGTGAACAATGCCTACCGGGACTGGCCAGCCTGGTGAAG
GTGTCAGAGCTGGCGCCTTCAAGTGGTGGTGCAGGAGGAGGGCGCAGTGGCCTCAGCCTCATCAAG
GAGACCTACCAGCTCCACAGGGACGACCCGAGGTGGTGGAGAAGCTGGGATGCTGCTGGTCCACCTG
GCTTCTATGAGGAGATCCTGCCGAGCTGGTGTCCAGTAGTATGAAGGCCCTGCTCCAGGAGATCAAG
GAGCGCTTACCTCCAGCCTGGTGTGAGTGACAGCAGCGCCTTACGAAACCAGGCCTCCCTCCAGGTGGA
AGCCCCAGCTGGGGTGCACCACGTCTGGGGACTGGAATAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map: □

ACCN: NM_153710

Insert Size: 2043 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

| | |
|-------------------------------|---|
| OTI Annotation: | This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_153710.4 |
| RefSeq Size: | 2794 bp |
| RefSeq ORF: | 2043 bp |
| Locus ID: | 169436 |
| UniProt ID: | Q8NE28 |
| Cytogenetics: | 9q34.2 |
| Protein Families: | Protein Kinase |
| MW: | 75.7 kDa |