

Product datasheet for **SC306664**

GSG1 (NM_153823) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GSG1 (NM_153823) Human Untagged Clone
Tag:	Tag Free
Symbol:	GSG1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC306664 representing NM_153823. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCAAGATGGAGCTCTCGAAGGCCTTCTGGCCAGCGGACACTCCTATCTGCCATCCTCAGCATG
CTACTACTCAGCTTCTCCACAACATCCCTGCTCAGCAACTACTGGTTTGTGGGCACACAGAAGGTGCC
AAGCCCTGTGCGAGAAAGGTCTGGCAGCCAAGTGCTTTGACATGCCAGTGTCCTGGATGGAGATACC
AACACATCCACCCAGGAGGTGTTACAATAACAACGGGAGACTGGGGATGACCGGTTCTCCTCCGGAGC
TTCCGGAGTGGCATGTGGCTATCCTGTGAGGAACTGTGGAAGAACCAGGGGAGAGGTGCCGAAGTTTC
ATTGAACCTACACCACCAAGAGAGAAATCCTATGGTTATCCCTGGGAACGCAGATCACCTACATC
GGACTTCAATTCATCAGCTTCTCCTGCTACTAACAGACTTGCTACTACTGGGAACCTGCCTGTGGG
CTCAAACCTGAGCGCCTTTGCTGCTGTTTCTCCTGTCTGTCTCAGGTCTCCTGGGGATGGTGGCCACATG
ATGTATTCACAAGTCTTCCAAGCGACTGTCAACTTGGGTCCAGAAGACTGGAGACCACATGTTTGGAA
TATGGCTGGGCCTTCTACATGGCCTGGCTCCTTACCTGCTGCATGGCGTCGGCTGTCACCACCTTC
AACACGTACACCAGGATGGTGTGGAGTTCAAGTGAAGCATAGTAAGAGCTTCAAGGAAAACCCGAA
TGCTACCACATCACCATCAGTGTTCCTCGGCGGCTGTCAAGTGCAGCCCCACCGTGGGTCTTTG
ACCAGCTACCACAGTATCATAATCAGCCATCCACTCTGTCTCTGAGGGAGTCGACTTCTACTCCGAG
CTGCGGAACAAGGATTTCAAAGAGGGCCAGCCAGGAGCTGAAAGAAGCAGTTAGTTCATCTGTAGAG
GAAGAGCAGTGTTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites:	Sgfl-MluI
Plasmid Map:	<input type="checkbox"/>
ACCN:	NM_153823
Insert Size:	981 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_153823.3
RefSeq Size:	2544 bp
RefSeq ORF:	981 bp
Locus ID:	83445
Cytogenetics:	12p13.1
Protein Families:	Transmembrane
MW:	36.7 kDa
Gene Summary:	May cause the redistribution of PAPOLB from the cytosol to the endoplasmic reticulum. [UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) uses an alternate splice site that causes a frameshift in the 3' coding region, compared to variant 1. The encoded isoform (2) has a distinct C-terminus and is longer than isoform 1. Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.