

Product datasheet for **SC306702**

TEDDM1 (NM_172000) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TEDDM1 (NM_172000) Human Untagged Clone
Tag:	Tag Free
Symbol:	TEDDM1
Synonyms:	EDDM9; Epdd1; HE9; HEL-S-45e; TMEM45C
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC306702 representing NM_172000. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
ATGATACTCAAAGGTTGTCTCCTGTATCCTTTGTGCTCTCCAGGAATAAGCAAAGATGTGCCAGGCTG
TGGAAATAGCCTATGGAGGATTACTAAAGATAGTACTGGCTCCCTCTTAACATTTTATGTGGTTCTC
TGTCTTGATGGCGGGATGGTGCTCATGAGGAAGCAAGTGCCATCAAGGTTTCATGTACCCAAAGAGTGG
CAGCACCTCACCATGTTTCATCCTCCTCACTCTTAATGGCTGTGTGGACTTCATGAGCAAGAATGTGCTG
CCTCAGAGGTGTGTGGGCTAGAAAAGGTACCCTGGTCTGATCATCTACGAGCTCCTGCTGCTGATG
GTGTCACATGTTAAAGATTGAGAAGGGTGGAGCTGCACGTTTATTCTCTGCTCATCTTGGTGGTGTTC
CTGCTGTTGCTGGTGTGACTGCAGAGCTGTGGGCTCCCAACATGTGTCATCTCCAGCTGATGGAGACC
TTTCTGATTCTGATGATGGGCTCCTGGCTGATGCAGGCAGGCTTTATTCTATACAGACCCGCTCTCTGGC
TACCCATGGCAGGACGATGACATCAGTGACATCATGTTTGTACCACTTCTTCTGCTGGCATGTGATG
ATCAATGCCTCATTCTGTTGGGAATCTATGGCTTCTCTTCTTTGGTATCATTGTTTCAGACCCAGC
TTGAAGCTGACTGGGCCAAAGAAGCTCCATATTATGCAAGCACTCCAGGACCCCTCTACAAGTTGCTA
CAGGAAGTGGAGCAGTCAGAGAAAGAGGACCAGGCTCTCTCTTCCAAAGAGCTCCCC**TGA**
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_172000
Insert Size:	822 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_172000.3</u>
RefSeq Size:	2500 bp
RefSeq ORF:	822 bp
Locus ID:	127670
UniProt ID:	<u>Q5T9Z0</u>
Cytogenetics:	1q25.3
Protein Families:	Transmembrane
MW:	31.3 kDa