

Product datasheet for **SC306725**

IL28A (IFNL2) (NM_172138) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	IL28A (IFNL2) (NM_172138) Human Untagged Clone
Tag:	Tag Free
Symbol:	IL28A
Synonyms:	IL-28A; IL28A
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL4</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for NM_172138 edited CCCTGGGTGACAGCCTCAGAGTGTTCCTCTGCTGACAAAGACCAGAGATCAGGAATGAA ACTAGACATGACTGGGGACTGCACGCCAGTGTGGTGCTGATGGCCGAGTGTGACCGT GACTGGAGCAGTTCCTGTGCGCCAGGCTCCACGGGGCTCTCCCGGATGCAAGGGGCTGCCA CATAGCCAGTTCAGTCCCTGTCTCCACAGGAGCTGCAGGCCTTTAAGAGGGCCAAAGA TGCCTTAGAAGAGTCGTTCTGCTGAAGGACTGCAGGTGCCACTCCCGCCTCTCCCCAG GACCTGGGACCTGAGGCAGCTGCAGGTGAGGGAGCGCCCATGGCTTTGGAGGCTGAGCT GGCCCTGACGCTGAAGGTTCTGGAGGCCACCGCTGACACTGACCCAGCCCTGGTGGACGT CTTGGACCAGCCCTTACACCCTGCACCATATCCTCTCCAGTTCGGGGCTGTATCCA GCCTCAGCCACGGCAGGGCCAGGACCCGGGGCCGCTCCACCATTGGCTGTACCGGCT CCAGGAGGCCCAAAAAAGGAGTCCCCTGGCTGCCTCGAGGCCTCTGTACCTTCAACCT CTTCCGCCTCTCACGCGAGACCTGAATTGTGTTGCCAGTGGGGACCTGTGTGTCTGACC CTCCACCAAGTCATGCAA
Restriction Sites:	Please inquire
ACCN:	NM_172138
Insert Size:	700 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	The ORF of this clone has been fully sequenced and found to be a perfect match to NM_172138.1.



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_172138.1</u> , <u>NP_742150.1</u>
RefSeq Size:	734 bp
RefSeq ORF:	603 bp
Locus ID:	282616
UniProt ID:	<u>Q8IZJ0</u>
Cytogenetics:	19q13.2
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane
Protein Pathways:	Cytokine-cytokine receptor interaction, Jak-STAT signaling pathway
Gene Summary:	This gene encodes a cytokine distantly related to type I interferons and the IL-10 family. This gene, interleukin 28B (IL28B), and interleukin 29 (IL29) are three closely related cytokine genes that form a cytokine gene cluster on a chromosomal region mapped to 19q13. Expression of the cytokines encoded by the three genes can be induced by viral infection. All three cytokines have been shown to interact with a heterodimeric class II cytokine receptor that consists of interleukin 10 receptor, beta (IL10RB) and interleukin 28 receptor, alpha (IL28RA). [provided by RefSeq, Jul 2008]