

Product datasheet for **SC306832**

SLCO6A1 (NM_173488) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLCO6A1 (NM_173488) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLCO6A1
Synonyms:	CT48; GST; OATP-I; OATP6A1; OATPY
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC306832 representing NM_173488.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTTTCGTAGGCGTCGCCCGCACTCTGGGAGCCAGGATGAAGTCTCAAGGGGAGTAGAGCCGCTGGAG
GCCGCGCGGGCCAGCCTGCTAAGGACAGGAGGCCAAGGGAACCCGAAGTCCTCGAAGCCCGGAAA
AAACACCGGTATCTGAGACTACTTCCAGAGGCCTTGATAAGGTTTCGGCGGTTTCCGAAAAAGGAAAAA
GCCAAGTCCTCAGTTTCCAAGAAGCCGGGAGAAGTGGATGACAGTTTGGAGCAGCCCTGTGGTTTGGGC
TGCTTAGTCAGCACCTGCTGTGAGTGTTGCAATAACATTCGCTGCTTCATGATTTTCTACTGCATCCTG
CTCATATGTCAAGGTGTGGTGTGGTCTTATAGATGTCAGCATTGGCGATTTTCAGAAGGAATATCAA
CTGAAAACCATTTGAGAAGTTGGCATTGAAAAGAGTTACGATATTTTCATCTGGCCTGGTAGCAATATTT
ATAGCATTCTATGGAGACAGAAAAAAGTAATATGGTTTGTAGCTTCTCCTTTTAAATAGGACTTGGGA
TCACTTTTATGTGCTTTTCCATCCATTAATGAAGAAAAAACAAGTAAGGTAGGAATTGAAGATATT
TGCGAAGAAAAAAGGTTGTCAGTGGTTCAGAGCAGTGGTATATCATTCCAATCAAAATACCTGTCT
TTCTTCATCCTTGGGCAGACTGTGACAGGAATAGCAGGAATGCCTCTTTATATCCTTGAATAACCTTT
ATTGATGAGAATGTTGCTACACACTCAGCTGGTATCTATTTAGGTATTGCAGAATGTACATCAATGATT
GGATATGCTCTGGGTATGTGCTAGGAGCACCAGTAAAGTCCCTGAGAATACTACTTCTGCAACA
AACACTACAGTCAATAATGGTAGTCCAGAATGGCTATGGACTTGGTGGATTAATTTTCTTTTGGCGCT
GTCGTTGCATGGTGTACATTAATACCATTGTCTGCTTTCCAACAATATGCCAGGTTCAACACGGATA
AAAGCTAGGAAACGTAACAGCTTCATTTTTTTGACAGCAGACTTAAAGATCTGAACTTGAACATAAT
ATCAAGGATTTATGTGCTGCTCTTTGGATTCTGATGAAGAATCCAGTGCTCATATGCCTAGCTCTGTCA
AAAGCTACAGAATATTTAGTTATTATTGGAGCTTCTGAATTTTTGCCTATATTTAGAAAATCAGTTT
ATATTAACACCCACTGTGGCACTACACTTGCAGGACTTGTTTTAATTCCAGGAGTGCACTTGGCCAG
CTTCTGGGAGGTGTCATTGTTTCCACATTAGAAATGTCTGTAAAGCCCTTATGAGATTTATAATGGTT
ACATCTGTGATATCACTTATACTGCTTGTGTTTATTATTTTGTACGCTGAATCCAGTGAATTTGCT
GGGATCAATGAAGATTATGATGGAACAGGGAAGTTGGGAAACCTCACGGCTCCTTGAATGAAAAATGT
AGATGCTCATCTTCAATTTATCTTCTATATGTGAAGAGATGATATTGAATATTTTCTCCCTGCTTT
GCAGGGGTGACATATTCTAAAGCACAAAACCAAAAAAGATGACTACAATTGTTCTTGCAATTAAGAA
GGATTAATAACTGCAGATGCAGAAGGTGATTTTATTGATGCCAGACCCGGGAAATGTGATGCAAAAGTGC
TATAAGTTACCTTTGTTCAATTGCTTTTATCTTTTCTACACTTATATTTTCTGGTTTTCTGGTGTACCA
ATCGTCTTGGCCATGACGCGGGTTGTACCTGACAACTGCGTTCTCTGGCCTTGGGTGTAAGCTATGTG
ATTTTGAGAATATTTGGGACTATTCCTGGACCATCAATCTTTAAATGTGAGGAGAACTTCTTGTATT
TTACGGGATGTTAATAAATGTGGACACACAGGACGTTGTTGGATATATAACAAGACAAAAATGGCTTTC
TTATTGGTAGGAATATGTTTTCTTTGAACTATGCACTATCATCTTCACTACTATTGCATTTTTCATA
TACAAACGTCGTCTAAATGAGAACACTGACTTCCAGATGTAAGTGTGAAGAATCCAAAAGTTAAGAAA
AAAGAAGAACTGACTTGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI
ACCN: NM_173488
Insert Size: 2160 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_173488.4</u>
RefSeq Size:	2689 bp
RefSeq ORF:	2160 bp
Locus ID:	133482
UniProt ID:	<u>Q86UG4</u>
Cytogenetics:	5q21.1
Protein Families:	Transmembrane
MW:	79.2 kDa