

Product datasheet for SC306848

C11orf42 (NM_173525) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C11orf42 (NM_173525) Human Untagged Clone
Tag:	Tag Free
Symbol:	C11orf42
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC306848 representing NM_173525. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
 ATGTTGGTGGGTACCCCAACCTGCTGACACTGGATGAAGCTGATGCCACCTGGACCCTCATCAAGGAT
 AAGGTCATCGAGGAGCACTTTGGGCCCAATGCAGTAGCAGTACCTTTCTGTGATGCAGCCTGCTAT
 GACCTACTGGGTGTGCTGGTAAACAGTCCCGCCAGCCCATACCCGCCTGGCTTTGCCAGGTCGGCAG
 GGCCGGAGGGCACTGAAACAGTGGGGCCACTACCAAGCCTCCTGGAGCAGGCAGGATCTGAGGGTGCC
 TTCGCCCACTGCACTCGGAATACTACCAAATGGCCGAGCAGAGAGAGCCTATGAAGAGACGCGAATG
 TTGGATGGACAGCCCTGCAAGATCCGCCTACATATGGGTGACCTGCGCAAGAAGTTGCCTTCTGTTG
 CTGCCACCAGGGCAGGTGAGCCTACAGCAGACTCTTCCCTGGCTCCGAAGCACCCACAGCATCTATGTC
 ATCTACCAGGTCTTCTTGTTCCTGGCTGCAGCTGGGGCTGACGTCTACAGCCCGTGAGCCCCAGCTC
 CTCCGGCTACTTCGGTCATTACCTGTTGCCTTCTCCTGCCTCAAGTTTTCACTGCAGTCTAAGGGCGTG
 CTGGGACCACAGAAGCCTCTCACTAAAGACCCATTGCCCATGGGGCCAACTGGGTGAGACCCAACCTC
 AGCATCATGCCGCTCTGGCCCCACATCAGCACCTGCTGATACAAGTGAAGCTGCTGATGTGCCCCCA
 CCTGTCCAGCCCCACCTACGCCACCTCCCAGGAAGGGCCAGAGGACAAACCCACCAGATTCTCCTAC
 AAGGGCCGAAACCCCTTCTGGAGGGGGCCCCAGATACTGTCAGAGAACTGGCTCTTCAGCCCCCGCAGC
 CCTCCACCAGGAGCCAGGGTGGGGGCCCGAGGACCCGACGGGCACTCCATGTCCTGCCCTGCTG
 CAGGGTCTATCCTCAGAGTTCGACAGTGACGACTGA
 ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_173525
Insert Size:	1002 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_173525.2</u>
RefSeq Size:	1149 bp
RefSeq ORF:	1002 bp
Locus ID:	160298
UniProt ID:	<u>Q8N5U0</u>
Cytogenetics:	11p15.4
MW:	36.4 kDa