

## Product datasheet for **SC306886**

### C18orf26 (DYNAP) (NM\_173629) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	C18orf26 (DYNAP) (NM_173629) Human Untagged Clone
Tag:	Tag Free
Symbol:	C18orf26
Synonyms:	C18orf26
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC306886 representing NM_173629. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGTTGCAGATATAAAGGGCAATGAACAAATTGAAAAATATTCTTGGAGAGAAGCTTGTGATACTGGC
AGCTCAAGAAATGGACAGAAAGCATGGAAAATACATATTGAACGTTGAGCACTCTGAAAACCGCCGCCA
ATCACACATCCAATGACCAAGAGGCTCACAGTCCATATGCTGGTGTCTACCTTCAAATGATATAACC
AGTGATGTCTCTCCAACCTAACTGGGGTCTGCGTGAACCCAGGAATCCTTGCACATTCAAGATGTCTA
CAGTCAGAATCCTGTAACACACAGGTAAGAATATTGCCCAATGACTGGTCTATGTGGAAAGTCTTC
CTGGCTTGTCTCTTAGCCTGTGTGATAATGACAGCAATTGGAGTACTTATAATATGCTTGGTGAATAAC
AAAGGATCGGCCAATTCCTCCATTGTTATCCAGCTATCCACAAATGATGGAGAGTGTGTGACTGTCAAA
CCTGGAACACCCTCTCCTGCTTGTCCACCTACAATGACCACCACTTCAACTGTACCTGCAAGTACAGCC
ACTGAATCTACAACCTCAACAGCTACAGCTGCCACCCTTCCACAGAACCTATAACTGTTGCACCTACC
GATCATTATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites:	Sgfl-MluI
ACCN:	NM_173629
Insert Size:	633 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_173629.2</a>
<b>RefSeq Size:</b>	2027 bp
<b>RefSeq ORF:</b>	633 bp
<b>Locus ID:</b>	284254
<b>UniProt ID:</b>	<a href="#">Q8N1N2</a>
<b>Cytogenetics:</b>	18q21.2
<b>Protein Families:</b>	Transmembrane
<b>MW:</b>	22.5 kDa
<b>Gene Summary:</b>	Plays a role in the regulation of cell proliferation. Promotes activation of the AKT1 signaling pathway. Promotes phosphorylation of AKT1 at 'Ser-473'. [UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.