

Product datasheet for **SC306933**

CXXC11 (RTP5) (NM_173821) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CXXC11 (RTP5) (NM_173821) Human Untagged Clone
Tag:	Tag Free
Symbol:	CXXC11
Synonyms:	C2orf85; CXXC11; Z3CXXC5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC306933 representing NM_173821.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGACCGGGCTGGGGCAGACATGTGGCCAGCACCTTCACCCTGGCCATGGCCGAGAGGAAGCCCCAG
GACGTCTGGGTTCTGCTACCTGAGCACAGCCTGGTCCCGGATGCCTGGACGGCGGTGGTGTCCAGTAC
CTGCTGGTGGGGCTCTCGAGGCTCCAGTGCAGTCACTGTCCGGGACCTGGGACTCGGCCCATGTGCAC
GTCTCTTCCACCTGTGGTGGGACAGGGCCAGCCACCGGGGGCTGGTGAAGATGCGCATCTGGGGCCAG
CGGTGCAGGCTGTGCCCCGACCCGGGACTGCCAGGTGAGGCCCGGGCGAGCAGCCCTTCTCAGC
AGGCTGGTCTTGACATCCTGCAGGACTGCTACGGGATGGCCCCGGCCAGCCGGCACCCAGGGAG
GCCTATGAGGGCTGTGTGAGGCTGTGAGCTGGGGTCTGCTTCTCCAGAAGGCCAGACCCGCC
TGGAGCGCAACGCCACAAAGGCAACTCCCCGCCACGGCTGGGTGGCACTGGCACCCTCTCCAGG
GGCAAACCGCTGTCCACCCTGGCGACGACCTTGGCAAGGGTGGCGTTGTCATCGCCATCCCCTTCTCC
CTTGTGGTACCAGCAATGACCAGGTGCCATCGCTGAGGGCCCTGCCCCCTCGGGGGCCTCTCTC
CCTGTGACTGGCAGCTGTGAGGCCTGGTCATCGGCCAGGGCTCCATCTTCTGTCTGGGGATTCAGTG
GCCATGCCTGGGGCAAAGGCTTCCCGGTGGCCATTGGAGACCCCTCTTCCACGGCCCCGGCTCCTC
GGCAGCAGCATCCAGACCTTCGAGCTCAAGGGCTTCTCTTCAAAGCGGGGCTCCCTCTGCAGCCG
GTTGGCGTGGCCAGGGCTGGGGCCCCATCTCCCTCAACAATGGCCTCGTCCCTGTGGGAAACACAG
CCAACCGTGTTCTACTGTGTGGGCTCTCGGCCAGCGGGAGGGCTCCCTCACCTTCCCCTCCTCCTC
ACCAGCATCTTACCAACACCCTCTCGGAGCCACCGATGGCCCTGTGGCACTAAAGAGGCTCCATC
ACCTTCCCCTTCTCTTACTGATGTCAAGGATGCCGTTGCTGAGGTGGTGAAGGCAACGGGAAGGAA
GGAGCGGCCAGGGCCTCGTCCAGTGGGTACGACGCCCTGCCAGAGACCAATGCTGGTGGCCTCCCC
TCCAGGTCAAGGGCTCCCTTGCCTCCCTTCCCTGCTGATGTCCAAGGCAAAGATGCCTTTACTGAC
ATCACTGAAGGCAAAGAGAAGGAAGGTGGCTGGTCACCGCGGGTCACGACGCCCTCTGGAGGCAAT
GCCGAGGGCCCCATCACGGTTAGTGAGGGCTGCATCACCATCCCCTTCGAGTCTTCGATGTCATAAAG
CGCAAGGGCGGTGGCCACGTTGCCTACGGCCCCAGGGCAATGGCTGCTTCTCCAAGGCTATTACCAG
AAGAGGCAGCTGAGGTCCAGGTTCCACAAGGCCGCTGTGGTGCCGCCGGGAGGAAGACGAGCGCCCT
GGCCGTGCTGCGTAGGCCGACGCCGAGCCCTACGAGGACTTCTGGATCTGGGTGTCATGACCGTG
TGCGTCTTCTGGCTGATGTGCATGTGTCGGCTGAACCCCGGGATCTACCCGACGAAAGTGA
ACGGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGCGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_173821

Insert Size: 1719 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_173821.2

RefSeq Size: 1798 bp

RefSeq ORF: 1719 bp

Locus ID: 285093

UniProt ID: Q14D33

Cytogenetics: 2q37.3

MW: 60.5 kDa