

Product datasheet for SC306952

BPIL3 (BPIFB6) (NM_174897) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BPIL3 (BPIFB6) (NM_174897) Human Untagged Clone
Tag:	Tag Free
Symbol:	BPIL3
Synonyms:	BPIL3; LPLUNC6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC306952 representing NM_174897. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
 ATGCTGCGGATCCTGTGCCTGGCACTCTGCAGCCTGCTGACTGGCACGCGAGCTGACCCTGGGGCACTG
 CTGCGGTTGGGCATGGACATCATGAACAGGTCCAGAGCGCCATGGATGAGAGTCATATCCTGGAGAAG
 ATGGCAGCCGAGGCAGGCAAGAAACAGCCAGGGATGAAACCTATCAAGGGCATCACCAATTTGAAGGTG
 AAGGATGTCCAGCTGCCCGTCATCACACTGAACCTTTGTACCTGGAGTGGGCATCTTCCAATGTGTGTCC
 ACAGGCATGACCGTCACTGGCAAGAGCTTCATGGGAGGGAACATGGAGATCATCGTGGCCCTGAACATC
 ACAGCCACCAACCGGCTTCTGCGGGATGAGGAGACAGGCTCCCGTGTTCAAGAGTGAGGGCTGTGAG
 GTCATCCTGGTCAATGTGAAGACTAACCTGCCTAGCAACATGCTCCCAAGATGGTCAACAAGTTCCTG
 GACAGCACCTGCACAAAGTCTCCCTGGGCTGATGTGTCCCGCCATCGATGCAGTCCTGGTGTATGTG
 AACAGGAAGTGGACCAACCTCAGTGACCCCATGCCTGTGGCCAGATGGGCACCGTCAAATATGTCTG
 ATGTCCGACACGACCCACAGCCAGCTACATCCAACCTGGACTTCAGTCTGTGGTGACAGCAAAAAG
 GGCAAAACCATCAAGCTTGCTGATGCCGGGAGGCCCTCACGTTCCCTGAGGGTTATGCCAAAGGCTCG
 TCGCAGCTGCTGCTCCAGCCACCTTCTCTCTGCAGAGCTTGCCTTCTGCAGAAGTCTTTTCATGTG
 AATATCCAGGATACAATGATTGGTGAGCTGCCCCACAAACCAAGACCCTGGCTCGCTTCATTCCT
 GAAGTGGCTGTAGCTTATCCCAAGTCAAAGCCCTTGACGACCCAGATCAAGATAAAGAAGCCTCCCAAG
 GTCATATGAAGACAGGCAAGAGCCTGCTGCACCTCCACAGCACCTGGAGATGTTTCGAGCTCGGTGG
 CGGAGCAAGGCTCCAATGTCCCTCTTTCTCCTAGAAGTGCACTTCAATCTGAAGGTCCAGTACTCAGTG
 CATGAGAACAGCTGCAGATGGCCACTTCTTTGGACAGATTACTGAGCTTGTCCCGAAGTCTCATCG
 ATTGGCAACTTCAATGAGAGGGAATTAAGTGGCTTCATCACCAGCTATCTCGAAGAAGCCTACATCCCA
 GTTGTCAATGATGTGCTTCAAGTGGGCTCCCACTCCCGACTTCTGGCCATGAATTACAACCTGGCT
 GAGCTGGACATAGTAGAATGCCCTGATGCTGGACTTGAAGCTGGG**TGA**
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTAAACGGCCGGC


[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_174897
Insert Size:	1362 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_174897.2</u>
RefSeq Size:	1362 bp
RefSeq ORF:	1362 bp
Locus ID:	128859
UniProt ID:	<u>Q8NFO5</u>
Cytogenetics:	20q11.21
Protein Families:	Druggable Genome, Secreted Protein
MW:	49.7 kDa