

Product datasheet for **SC306991**

DDX51 (NM_175066) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DDX51 (NM_175066) Human Untagged Clone
Tag:	Tag Free
Symbol:	DDX51
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF:	<p>>NCBI ORF sequence for NM_175066, the custom clone sequence may differ by one or more nucleotides</p> <pre> ATGGCGCTGTTCTACGTCGCGCGGTACCCGGGCCCGATGCGGCAGCTGCGGCGGGGCCG GAGGGCGCGGAGGCCGGGGCGCACGGCAGGGCCCGCGCTGCTCGAGCGGTGCAGAGC CGGGCCCGCAACGGCAGCAGCAGCGGGAGCCCGCGCAGACCAGGCGGTGCATCGACC GAGCCGGCGACCAGGAGGCGACGGCGCCCGGGCGGGCGGGTGAACGACGCGGAG CCGGGAGCCCGGAGGCGCCGAGGAAAGCGACGGAAGGCGGACGGCGAGGACGCGGGC GCAGAAAGCAACGAGGAGGCGCCAGGGGAGCCAGCGCAGGGAGCAGCGAGGAGGCGCCA GGGAGCCAGCGCAGGGAGCAGCGAGGAGGCGCGGGGAGCGCAGCACCAGCGCCAGC GCCGAGGCGGCCAGATGGACCGGCCCTGGAGGAGGCGCCGGACCCCTGGTCCCCGGC CTGGTGTGGGGGGTTCGGGAAGAGGAAGGCGCCGAAGTCCAGCCTTCTGCCAAGG TGGCTGGTGAGCCTAAGTGTGTCAGAAGGAATGTACCGAAGACCTGGTTCCTATCGAG GACATCCCTGACGTCCATCCTGACCTGCAGAAGCAGCTGCGGGCACACGGCATCTCGTCC TACTTTCCAGTCCAGGCAGCTGTGATTCTGCCCTCTGGAGAGCGCAGCCTGTGGGTTT CTGGTGGGCGAGGTTGGCTACCGCCTAGCGACCTCTGTGTTTCTGCCCCAACAGGCAGT GGGAAGACTGCGCTTGTGCATCCCTGTGGTGCAGGCCCTGCTTTCGAGAGTGGTCTGC CACATCCGTGCCCTGGTTGTGCTGCCACCAAGGAGCTGGCCAGCAGGTGAGCAAAGTT TTCAACATCTACACAGATGCCACACCTCTGAGAGTCTCCCTGGTACGGGACAGAAGTCT CTGGCCAGGAGCAGGAGAGCCTCGTCCAGAAAACAGCTGATGGGTACCGCTGCTTGGCT GACATCGTGGTAGCCACCCCGGCCCTGGTGGACCACATCGACCAGACCCAGGATTC AGCCTCCAGCAGTCCGCTTCTGATTATCGACGAGGCTGACCGGATGATTGACAGCATG CATCAGTCTGGTGCCTGGGTGGTGGCGCCGCTTCCAGAGCGAGGACCCCGCGGAC CCCTGTGCCCTGCTCCAGCGAAGCAGGCCAGGCTGTGACAGCCGCCAGCACCTGCTGT CCCCAGATGCCCTGCAGAAGCTGCTTCTCAGCTACTCTGACCCAGAACCCTGAAAAG CTGCAGCAGTGGCCTCCACCAGCCCGGCTTTCTCCACAGGCTAGCACACAGGGGC CTGGAAGATACAGATGGGGACGGGGATTGGGGAAGTATGCCTTCTGTTGGGCTCACG CACCCTACGTGCCCTGCAGCCTCAGCTTAAGCCGCTGGTGTGCTGACCTGGTCTG GAGATGGGCTTCTCGAGGTTCTGCTTCACTAACTCCCGAGAGAACTCCACAGGCTC TTCCTGCTGGTGAAGCTTTGGGGGTGTGGACGTGGCTGAGTTCTCCTCGCGTACGGG CCTGGCCAGAGGAGGATGATCCTGAAGCAGTTTGAACAGGGGAAGATCCAGCTGCTCATC AGCACGGACGCCACCGCGGAGGCATCGACGTGCAGGGTGTGGAGCTGGTGGTGAACAC GACGCCCCAGTACCTGAGAACCTACGTGCACCGGGTTGGGAGGACAGCTCGCGCTGGG AAAAGTGGACAGGCTTCACTGCTCCTGAAAGTGCAGGAGAGGAGATTCTCCGAATG CTAACTGAAGCTGGGGCACCTGAGTTGCAGCGGCACGAGCTCTCCAGCAAGCTGCTGCAG CCGCTGGTTCCTCGGTACGAGGAGGCCCTGTCCAGCTGGAGGAGTCTGTCAAGGAAGAG CGCAAGCAGAGGGCGGCCTAG </pre>
Restriction Sites:	Please inquire
ACCN:	NM_175066
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_175066.2](#), [NP_778236.1](#)

RefSeq Size: 2802 bp

RefSeq ORF: 2001 bp

Locus ID: 317781

UniProt ID: [Q8N8A6](#)

Cytogenetics: 12q24.33

Gene Summary: ATP-binding RNA helicase involved in the biogenesis of 60S ribosomal subunits.
[UniProtKB/Swiss-Prot Function]