

## Product datasheet for SC307129

### FADS6 (NM\_178128) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	FADS6 (NM_178128) Human Untagged Clone
Tag:	Tag Free
Symbol:	FADS6
Synonyms:	FP18279
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC307129 representing NM_178128. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAACCCACGGAGCCGATGGAACCTACGGAGCCCATGGAACCTACGGAGCCCATGGAACCTACGGAG
CCCATGGAACCTACGGAGCCCATGGAACCCGGCGGAGCGCGCACCGTGGGGGCGAGGCGCTGCTGCGG
GAGCTGGAGGTGCTGGTGCAGGACGTGGTGAAGACGAGCTCCTGGTGGGAGCGCCACGGCGTGGACTGC
GCCATCCTCGCGCTCAGCCTCTTCGCTTGCCTTGCCTGCGGCGAGGCTTCTGTGCCTGCGCTGGGAGAATGCCCTG
GTCTTTGCATCCGGCATCACCATCTTGGGTGTGTGCCACTACACACTCACTGTCAAGGGCAGCCACCTG
GCCACTCATGGGGCCCTCACCGAGTCAAACGCTGGAGCAAGATCTGGCTGCTTTTCTTTGTGGAGGTG
TGCACAGCCTTCACTGCAGAGCACGCCACGCATGGGCACGTCAAGATGCACCATGCCTACACCAAGGTG
GTGGGCTGGGGGACTCCAGCACGTGGAGGCTGCCTTGCCTCAACCGCTATGTCTACATGTTCTTGTCT
CCTTTCTCCTCCCATCGCCACTCCACTGGTGGCTGTCGAGCGGCTGAGGAAGGTGGAGCTCGGGACA
GCCCTGCGGACGCTGGCCCTGATTTCTCTGGGCCCTTATTCTCACTACTGGCTGCCTGAACGTGTCA
GGCTTCAAGAACCCAGCTCAGCCCTGGGCTGCATGTTCTCACCAGATCCCTGTTGGCCACCCCTAC
CTCCACGTCAACATCTTCCAGCACATCGGACTGCCATGTTCTCCCGGACAACAAGCCCGTCGGATT
CACATGATGAGCCTGGGGTGTAACTGGCCCGGCTGCCGTGCTGGACTGGGCGTTCGGCCACTCG
ATCATCAGCTGCCATGTGGAACACCATCTATCCAGGCTCTGATAACATGTGCCTGAAGGTGAAG
CCGTGGTGTCCAGTTCCTACGTGAGAAGCAGCTACCGTACAACGAGGACTCATACCTGGCTCGCTTC
CAGCTGTTTCTCGTGCCTATGAGGAATTCATGGTGCAGGCCCCACCCATCACTGAGCTTGTGGGCTG
TAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

<b>Plasmid Map:</b>	□
<b>ACCN:</b>	NM_178128
<b>Insert Size:</b>	1107 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_178128.5</a></u>
<b>RefSeq Size:</b>	2352 bp
<b>RefSeq ORF:</b>	1107 bp
<b>Locus ID:</b>	283985
<b>UniProt ID:</b>	<u><a href="#">Q8N9I5</a></u>
<b>Cytogenetics:</b>	17q25.1
<b>Protein Families:</b>	Transmembrane
<b>MW:</b>	41.8 kDa