

Product datasheet for **SC307145**

PILRB (NM_178238) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PILRB (NM_178238) Human Untagged Clone
Tag:	Tag Free
Symbol:	PILRB
Synonyms:	FDFACT1; FDFACT2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC307145 representing NM_178238. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGGTCGGCCCTGCTGCTGCCCTGCTGCTCCTGCTGCAGCCGCCAGCATTCTGCAGCCTGGTGGC
TCCACAGGATCTGGTCCAAGCTACCTTTATGGGGTCACTCAACAAAACACCTCTCAGCCTCCATGGGT
GGCTCTGTGAAATCCCCTTCTCCTTCTATTACCCCTGGGAGTTAGCCATAGTTCCCAACGTGAGAATA
TCCTGGAGACGGGGCCACTTCCACGGGCAGTCCTTCTACAGCACAAGGCCGCCTCCATTACAAGGAT
TATGTGAACCGGCTCTTCTGAAGTGGACAGAGGGTCAGGAGAGCGGCTTCTCAGGATCTCAAACCTG
CGGAAGGAGGACCAGTCTGTGTATTTCTGCCGAGTCGAGCTGGACACCCGGAGATCAGGGAGGCAGCAG
TTGCAGTCCATCAAGGGGACAACTCACCATCACCCAGGCTGTACAACCACCACCACCTGGAGGCC
AGCAGCACAACCACCATAGCCGGCCTCAGGGTCACAGAAAGCAAAGGGCACTCAGAATCATGGCACCTA
AGTCTGGACACTGCCATCAGGGTTGCATTGGCTGTGCTGCTCAAACTGTCATTTTGGGACTGCTG
TGCCCTCCTCCTGTGGTGGAGGAGAAGGAAAGGTAGCAGGGCGCCAAGCAGTGACTTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites:	Sgfl-Mlul
Plasmid Map:	<input type="checkbox"/>
ACCN:	NM_178238
Insert Size:	684 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_178238.3
RefSeq Size:	1438 bp
RefSeq ORF:	684 bp
Locus ID:	29990
UniProt ID:	Q9UKJ0
Cytogenetics:	7q22.1
Protein Families:	Druggable Genome, Transmembrane
MW:	25.5 kDa
Gene Summary:	The paired immunoglobulin-like type 2 receptors consist of highly related activating and inhibitory receptors that are involved in the regulation of many aspects of the immune system. The paired immunoglobulin-like receptor genes are located in a tandem head-to-tail orientation on chromosome 7. This gene encodes the activating member of the receptor pair and contains a truncated cytoplasmic tail relative to its inhibitory counterpart (PILRA), that has a long cytoplasmic tail with immunoreceptor tyrosine-based inhibitory (ITIM) motifs. This gene is thought to have arisen from a duplication of the inhibitory PILRA gene and evolved to acquire its activating function. [provided by RefSeq, Jun 2013]