

Product datasheet for **SC307165**

LCE1C (NM_178351) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LCE1C (NM_178351) Human Untagged Clone
Tag:	Tag Free
Symbol:	LCE1C
Synonyms:	LEP3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC307165 representing NM_178351. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
 ATGTCCTGCCAGCAGAGCCAGCAGTCAGTCCAGCCCCCTCCCAAGTGCACCCCAAGTGCCTCCCAAG
 TGCCCCACCCAAAGTGTCCCCAAAGTGTCCCCCTAAGTGCCCTCCTGTCTTCTCCTGTCAGTGTG
 AGCTCCGGAGGCTGCTGTGGCTCCAGCTCTGGGGCAGCTGTGGCTCCAGCTCTGGGGGATGCTGCAGT
 TCTGGGGAGGTGGCTGCTGCCTGAGCCACCACAGGCGCCGTAGGTCCACTGCCACAGACCCAGAGC
 TCTGGCTGCTGCAGCCAGCCCTCGGGGGCTCCAGCTGCTGTGGCGGGGAGTGCCAGCACTCTGGA
 GGCTGCTGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_178351
Insert Size:	357 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.


[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_178351.3</u>
RefSeq Size:	695 bp
RefSeq ORF:	357 bp
Locus ID:	353133
UniProt ID:	<u>Q5T751</u>
Cytogenetics:	1q21.3
MW:	11.5 kDa
Gene Summary:	<p>Precursors of the cornified envelope of the stratum corneum.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1). Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.</p>