

Product datasheet for SC307168

LCE1F (NM_178354) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LCE1F (NM_178354) Human Untagged Clone
Tag:	Tag Free
Symbol:	LCE1F
Synonyms:	LEP6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC307168 representing NM_178354. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCTTGCCAGCAGAGCCAGCAGCAGTCCAGCCCTCCCAAGTGCCTCCCAAGTGCCTCCCAAG
TGCCCCACACCGAAGTGCCTCCCAAGTGTCCCTTAAGTGCCTCCTGTCTCTTCTGCTGCAGCGTC
AGCTCCGGAGGCTGTGTGGCTCCAGCTCTGGGGCTGCTGCAGCTCTGGGGTGGTGGCTGCTGCAGC
TCTGGGGAGGCGGCTGTTGCCTGAGCCACCACAGACGGCGTAGGTCCCACCGCCACAGACCCAGAGC
TCTGACTGCTGCAGCCAGCCCTCAGCGGCTCCAGCTGCTGCGGAGGGGGCAGTGCCAGCACTCTGGA
GGCTGCTGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:

ACCN: NM_178354

Insert Size: 357 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_178354.2
RefSeq Size:	598 bp
RefSeq ORF:	357 bp
Locus ID:	353137
UniProt ID:	Q5T754
Cytogenetics:	1q21.3
MW:	11.7 kDa
Gene Summary:	Precursors of the cornified envelope of the stratum corneum.[UniProtKB/Swiss-Prot Function]