

Product datasheet for **SC307187**

FAM83C (NM_178468) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FAM83C (NM_178468) Human Untagged Clone
Tag:	Tag Free
Symbol:	FAM83C
Synonyms:	C20orf128; dj614O4.7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC307187 representing NM_178468.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGATCGCC
ATGTTCCGAGGCCCGGGGCTGGGTCTGGGAGCCAGGGCATGGCGGGACCCCTGCGGGCGGGGTG
GAAGAGCTGAAGCTGCCGTGGTGGCGGGAGAGCTACCGCTGGTCTGCGGCACAGCGAGGCGGCTCGG
CTGGCGGCCGACGCCCTCCTGGAGCGGGTGAGGCTGCCACCTGCGGGTCATCTCCGAGGAGCGGGAG
CTGCCCTTCTGAGCGCCCTGGATGTGGACTACATGACCAGCCATGTGCGCGGGGGCCCTGAGCTCAGC
GAGGCTCAGGGGAGGAGGCTCCGGGCCAGACCGCTCAGCCTGCTCTGAAGTCACCTCAGGGACT
TACTTCCCCATGGCCTCTGACATAGACCCCCAGACCTGGACCTGGGCTGGCCGAGGTGCCACAGGCC
ACAGGCTTCAGCCCCACCCAGGCTGTGGTCCACTTCCAGAGGGACAAGGCCAAGAACATCAAGGACCTG
CTGCGCTTCTTTTCAGCCAGGCCACACGGTGGTGGCTGTGGTATGGACATATCACTGACATGGAG
CTTCTGTGTGACCTCATGGAGGCTCAAGCCGGCGTGGTGTCCCTGTGTACCTGCTCCTGCCAGGAG
CACCTGAGGCACTTCTGGAGATGTGCTACAAGATGGACCTCAATGGGGAGCACCTGCCGAACATGCGT
GTGCGGAGCACGTGTGGGACACATACTGACGCAAGGCTGGCCGCGCTTCACGGGGCAGGCCCTGGAG
AAGTTCGTCTCATTGACTGTGAGCAAGTGGTGGCGGGCAGTTACAGCTTCACTGGCTTTGACGCCAG
GCCCACTAGCATGGTGTGACGTGAGGGGCCGATCGTGGAAGACTTTGACCGGAGTTCGGCTGT
CTGTACGTGAGTCGACGCTGTGGAGGCTTCTGTGGCGGTGAGGACCGCTGTCTCCCGGGCACTG
CGTCTCCCCCTGTGGCCCTAGCCTTACGGCTGTATGCCAAGCCCCACGTGTCCTGCCCTCCAGC
ACCAGCCTCAGCAGCATCAAGCAGTCACCGCTTATGGTGCCTCCTACCTCGCTCTACCAGGAGGT
GGTGATTGCAGTGATACGGGTGTGGTGTCTCGTCCCTGGTCTGCGCGCGGTGAGGCCAGTGGCCAG
CCCTCCCTACATCGCAACTGTGACACCTAACCACGGCTCCCTCCTGGGCTCTATAGGGCCAATCTC
GGCAAGCTAGGGGCATACCCATGGTCCAGTCTCCCTGCCCTCAACCATAATAGTACCAGCCCCCTTA
ACCTTGGCAGTGGGGTCACCTCTGCTTCTCGCTCCCGGCCCTCCTCCAGTTCCATCGGGGTGCCCA
GCTCTGTCCCGTTCCAGAGAATGGGCTCCAGGAAGCCAAGAGCCAGCCCCCTGCGGGTGCATGG
GTACCTGGCACAACCCTGGAGACAGTGGAGGAGAAGGAGAAGAAGGCATCTCAAGTCAGAGCCGTGGC
CAGCTGGATCTCCTGTCCCTTCCCCAGAGCCGAGAAGTGGGAGACCCTGACTCTGGGGTTACCCCC
AACTCAGGCCCTTCCGGCTGGCAGCAGGCCCCAGAGGACAGGAGGTGTCCCAAGCCAGGCCGAC
AGCCAGCTGGATCTCCTGTCCCGAGCCCTGGGTACTGGGGTGCCCTGAGTTGGGTTCCTCAGACCT
GGTGATCGGGCCCTGGAGGACAGGAGGCTGTCCCTAAACCAAGCCGTGGCCAATCAGACCTCCTGATG
CAGTACCCCAAGGCCAGGGTTCAGAGTGCCCTTGAAACCAACTCCTCAGCCAGACCTGCCAGACGG
GCACCAGATGAGCGGCGGACAGCCCTGGGGCACAGCCAGCTGGACCTCATCACAAGTTCGCGCCATTC
CGTGGTGAGGGGCTGGGCCAATGGTCTCCCGATATCAAGCCCTGCTCGCACGGCTGGAGCTGGGTCT
GGGGATGAGAAACGGCTAACCTGGGCCACAGCAAGCTGGACCTCATCACAAGTATCATCAGTTGCAC
GGGGCAGGCAGGGAAGTGAAGCTGGGGTCCCAAGGTGGCCATCTCAATGGTGGTAACAGTGACCTG
GTCAGGGATGAGAAACGGCTGACCCTGGGTACAGCAAACTGGACCTCATCACTAAGTACAACAAGTCC
AAGTTCAAGCAGCTCCGAAGCCGCTTTGAGTCCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI
ACCN: NM_178468
Insert Size: 2244 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_178468.5</u>
RefSeq Size:	3166 bp
RefSeq ORF:	2244 bp
Locus ID:	128876
UniProt ID:	<u>Q9BQN1</u>
Cytogenetics:	20q11.22
MW:	81.1 kDa
Gene Summary:	This gene encodes a member of the family with sequence similarity 83 protein family. The encoded protein may be involved in regulating MAPK signaling in cancer cells. [provided by RefSeq, Mar 2017]