

Product datasheet for SC307219

TBC1D26 (NM_178571) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: TBC1D26 (NM_178571) Human Untagged Clone
Tag: Tag Free
Symbol: TBC1D26
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >SC307219 representing NM_178571.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAGATGGATGGGGACCCGTATAACCTGCCTGCCAGGGGCAAGGCAATATCATCATTACTAAGTAT
GAGCAGGGACACCGAGCTGGGGCAGCAGTGGACTTGGGGCATGAGCAGTTGATGTCAGAAAATACACC
AATAACCTCGGGATTGTGCATGAGATGGAGCTGCCCGCGTCAGTGCCCTGGAGGTGAAGCAAAGACGC
AAGGAAAGTAAACGTACCAACAAGTGGCAAAAGATGCTTGCACTGGACAAAATATAGGAGCACCAAG
AAGCTGTCTCAAAGAGTATACAAAGTCATTCCCCTGGCGGTACGGGGCCGGCGTGGTCACTTTTGCTA
GATATTGACAGAATCAAGTCCCAGAACCCAGGCAATATAAGGTCATGAAGGAGAAGGGCAAGAGGTCC
TCCAGAATCATCCACTGCATCCAGCTAGATGTCAGCCACACCCTGCAGAAACACATGATGTTTATACAA
AGATTGCGAGTCAAGCAGCAGGAATTATGTGACATCCTCGTGGCCTATTCTGCATATAACCTGAGGTG
GGCTACCACAGGGACCTGAGCCGCATCACTGCCATCCTCCTGTGTCTGCCAGAGGAAGATGCTTTC
TGGGCGCTTACCCAGTTGCTCGCTGGTGAGAGGCACTCCCTGTGGTATTCTACAGCCAAATACTGCCT
GGCTCGAGAGGCTCCTATCGCACCAGGAGCAGGTGCTGCACAAGTCCTCCCAAAGATCATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_178571
Insert Size: 753 bp
OTI Disclaimer:

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_178571.4</u>
RefSeq Size:	1883 bp
RefSeq ORF:	753 bp
Locus ID:	353149
UniProt ID:	<u>Q86UD7</u>
Cytogenetics:	17p12
MW:	28.9 kDa
Gene Summary:	May act as a GTPase-activating protein for Rab family protein(s).[UniProtKB/Swiss-Prot Function]