

Product datasheet for **SC307235**

OST beta (SLC51B) (NM_178859) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	OST beta (SLC51B) (NM_178859) Human Untagged Clone
Tag:	Tag Free
Symbol:	OST beta
Synonyms:	OSTB; OSTBETA
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_178859, the custom clone sequence may differ by one or more nucleotides

```
ATGGAGCACAGTGAGGGGCTCCCGGAGACCCAGCCGGTACTGTGGTACCCAGGAGCTGCTGGAAGAGA  
TGCTTTGGTTTTTTCGTGTGGAAGATGCATCTCCCTGGAATCATTCCATCCTTGCCCTGGCAGCTGTGGT  
GGTCATTATAAGCATGGTCCTCCTGGGAAGAAGCATCCAGGCAAGCAGAAAAGAAAAGATGCAGCCACCA  
GAAAAAGAACTCCAGAAGTCTGCATTTGGATGAGGCCAAGGATCACACAGCCTAAACAACCTAAGAG  
AACTTTGCTCTCAGAAAAGCCAACTTGGCCAGGTGGAACCTTGAGTTAAAAGAGAGAGATGTGCTGTC  
AGTTTTCTTCCGGATGTACCAGAACTGAGAGCTAG
```

Restriction Sites:	Please inquire
ACCN:	NM_178859



[View online »](#)

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_178859.2](#), [NP_849190.1](#)

RefSeq Size: 454 bp

RefSeq ORF: 387 bp

Locus ID: 123264

UniProt ID: [Q86UW2](#)

Cytogenetics: 15q22.31

Protein Families: Transmembrane

Gene Summary: Essential component of the Ost-alpha/Ost-beta complex, a heterodimer that acts as the intestinal basolateral transporter responsible for bile acid export from enterocytes into portal blood. Efficiently transports the major species of bile acids. Modulates SLC51A glycosylation, membrane trafficking and stability activities.[UniProtKB/Swiss-Prot Function]