

Product datasheet for **SC307418**

ZSWIM2 (NM_182521) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZSWIM2 (NM_182521) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZSWIM2
Synonyms:	MEX; ZZZ2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC307418 representing NM_182521.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTTCGCCGAGGCTATAAGGCCTCTGAAAGGCGAAGACACTTGAGCGAGAGGCTCAGCTGGCACCAA
GACCAGGCGCTGAGTAGCAGCATCTACCTCCTACGAGAGATGGGCCCACTGGCTTCCTGCTGAGGGAG
GAGGAGCCGGAATACATGGATTTCCGAGTTTTCTAGGAAATCCTCACGTTTGTAACTGTTCCACATTT
CCGAAAGGAGGGGAACCTTTGTAAGCATATCTGCTGGGTCTTGTGAAAAAATTCAAGCTTCCAAGGAAC
CATGAATCTGCTTTACAACCTGGGTCTTGAGAAAGAGAGATAAGTGACTTGCTTCGGGGGATACATCGA
GTTCAAACCTCCCAACCAGGAACAAATGACGAAAAATGAACATGTTGAAGAAGATGGGTACATTAACAG
AAGGAAATTGATTCAGAGGATATCTGCTCTATTTGTCAAGAGCTACTTTTAGAGAAAAAGCTTCCTGTC
ACCTTTTGCAGGTTTGGCTGTGGCAATAGTATTCATATAAAATGCATGAAGATCTTAGCTAATTATCAG
AGTACATCAAACACTTCCATGTTGAAATGCTCTGTGCAGGAAAGAGTTTGCACCATTAACACTGATT
TTGGAGGAATTCAAAACTCTAGCAAACTAGTAGCTGCAGCAGAAAAAGAGAGACTGGACAAACACCTT
GGGATTCCTGTAAATACTGCAACAGTTTCCAATTGAGGGGAAGTGTATAAGTGATACCGAATGCATA
GAATATCACTTATGCCAGGAATGTTTTGATAGCTGCTGCCATCTTTCCACACGTTTACATTTCTGTGAG
AAAAGAAACAAAAATGGAGATCACTAGAAAAAGAGCAGATGAAGTTGTAATAACATAGATACTAAA
AATGAGATTGAAGAAAAGATGTCACATTTTCAAGAAAAGCAAGGCCAAGTTTACACACCAAAACACATT
GTAAGATCACTGCCTCTCCAACCTGATTACTAAGAATAGTAAGCTGCTTGTCCAGGCTACCAAGTGTCTA
CTTTGTTTGAAGGCATTTTCATCTTGGTCAACATACAAGATTGCTACCATGTAATCACAAGTTTACAGG
AAGTGTATTGACAACCTGGTTATTCACAAGTGCAATTCATGCCCTATTGATGGACAAGTTATATATAAC
CCTTTAACTTGGAAAAATTCAGCAGTGAATGGACAAGCACATCAGTCTGTTTCAAACAGAGACATCATT
CATCTATCAAAGCAGAAAGAACCAGATCTTTTTATTCTGGTACTGGATTAGTCTTAAACAAAAATAGA
CTTGAATTTTACCTAGCATACCTCAGTGAATTTTGTGAATTGAATACACCTCAAAGCCCAAAAGAT
GCCTATGAAAATACAACAATAGATAATCTATGCTCTATCAAATTAGATAATTCAAATTCAAAAAATTA
ACCTATGATTATAAAATTAGCCAACATTTTCCAGGTATCTTCAAGATTACCCACTGTGTCATTTGGG
AAAATACCATCTCAAACACTGCTTCCTCTATTGTTTATAAGAATATTGTGTGCCACTGCAATGGAA
AGTCCATGCATCAGTGGAAAAATTCACACTAGTCTAAGCCGGATGACCAAGGCTGTAATGTAATAAC
CACAACCTAAAGAAGACTCCTGCCACTAAAATAAGAGAGGACAACAAGAGATCACTTTACTTCCAGAG
GATTTCAATCTTATTGTCAATTGGAGCACAGCTAACTTAGTTTGTCTAAAAGGTATAGTAAGTGTATG
GGGAAATTACACGAAAATGTAGTCATCTATCAAGACAGCCTGTGTCTCACTCTGTAATACAAAAAGT
ACTGAGCTATCTTTAATAATAGAAGGAGTTCAATTGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_182521

Insert Size: 1902 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_182521.2</u>
RefSeq Size:	2446 bp
RefSeq ORF:	1902 bp
Locus ID:	151112
UniProt ID:	<u>Q8NEG5</u>
Cytogenetics:	2q32.1
MW:	72.7 kDa
Gene Summary:	E3 ubiquitin-protein ligase involved in the regulation of Fas-, DR3- and DR4-mediated apoptosis. Functions in conjunction with the UBE2D1, UBE2D3 and UBE2E1 E2 ubiquitin-conjugating enzymes.[UniProtKB/Swiss-Prot Function]