

Product datasheet for **SC307467**

CCDC37 (CFAP100) (NM_182628) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CCDC37 (CFAP100) (NM_182628) Human Untagged Clone
Tag:	Tag Free
Symbol:	CCDC37
Synonyms:	CCDC37; MIA1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC307467 representing NM_182628.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGATCGCC
ATGTCTGAGATACCGTCCACTATAGTCTCCAAGAACATGACCAATGACAAGAACAGCCTGGAATCCATG
AACATCAGCTCTTCTTCAAGCACTGAAGAGAACCACAAAGAAACAGGCAAGAAAAACGAAGAACATGGT
CCTGACCCTTCAGCGAACCCTTTCCACTTATCTGGGGATGTGGATTTCTTCTTGCTCAGAGATCAGGAG
CGGAATAAGGCTCTCTCCGAACGGCAGCAGCAGAAGACGATGCGGGTGCACCAGAAGATGACCTACTCC
TCGAAAGTGTCGGCTAAGCACACCAGCCTGCGCGGCAGCTGCAGCTGGAGGACAAGCAGGAGGACCTG
GAGGCGCGCGCCGAGGCCGAGCATCAGCGCGCCTTCCGCGACTACACGACCTGGAAGCTCACCTTGACC
AAAGAAAAGAATGTGGAGCCTGAGAACATGAGTGGCTACATTAAGCAGAAGCGGCAATGTTCTCTCTC
CAGTATGCCCTGGATGTCAAGCGGAGAGAGATCCAGCGGCTGGAGACGCTGGCGACCAAGAGGAGGCC
AGGCTGGAGCGGGCCGAGAAATCCCTGGAGAAGGACGCCGCTTGTTCGACGAGTTCGTACGGGAGAAT
GACTGCAGCTCCGTGCAGGCCATGAGAGCGGCTGAGAAGGAGACCAAGCAAGATAGAGAAGATCCTT
GAGATCCGGGACCTCACCACCCAGATTGTTAATATCAAAAGTGAGATCTCCAGATTTGAAGACACTCTG
AAGCATTACAAGGTCTATAAGGATTTCTATACAAGCTGTCGCCCAAGGAGTGGCTTGAAGAACAGGAA
AAGAAACACTCGTTTCTCAAAAAGGCCAAGGAGGTCTCCGAGGCTTCCAAAGAGAGCAGTGTTAACTCC
ACACCAGGGGACAAAGGACCAGGGATCAAGGGCAAGGCGAGCTCCATGTGGGCCAAAGAGGGTCAGGGT
ACAAAGAAGCCCTGGAGGTTTCTGCAGACGATGCGGCTGGGGCGGAGCCCGTCTTACCTGAGCAGCCCC
CAGCAAGGCAGCCAGCCAGCGAGTCCAGCGGTGGCGACTCCAGAGGGTCGAACCTCTCCATCCCCCCC
ACGCAGGAGGACACCGACAGCGATGGGGAGGAACACAGCTGTACTTCACGGAGCCCCAGCAGCTCCTG
GATGTCTTCCGAGAGCTGGAGGAGCAGAACCTGTCGCTGATCCAGAACAGCCAGGAGACGGAGAAGACC
CTGGAGGAGCTGAGCCACACCCTGAAACACACCAGATCCGCATGGACAGGGAGGTCAACCAGCTGAAG
CAGTGGGTCAACCAATGATGATGTCCATCACCAAGGAGGAGGACACAGCAGCTGAGCTGGAGCTCAAA
GCCCAGTCTTCCACTTCGGCGAGTACAAGGGCGATCAGCAGGATAAGCTGCTAGAGAGCCTGAAGTGC
AAGGTGCTGGATGTGTACCGGCACTGCACCGGCACCCAGCAGGAGGCCAACCTGGGCACCGTGCAGATG
CTGACCATCATTGAGCACCAGCTGGATGAGCTGCTAGAGAACCTGGAGCACGTGCCCCAGGTCAAGATC
GAGCAGGCCGAGAGGGCAAAGGAGAAGGAGCGGCGCATCAGACTTCGAGAAGAGAAGCTCCAGATGCAA
AAGATCTACAGGAGGAGCATCTGCAGCGGGCCCGGGCGCGCCAGGCTGAGATCAAGAAGAAGAGA
GGCAGGACACTGGTATGCCGCTCAGACCCCCAGCCACAGGATCAAACAACAGTCTGAGCACACACTG
ATGGACAAGGAGGAGGAGGAGCTGCTATTTTCTTTACTTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_182628

Insert Size: 1836 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_182628.2</u>
RefSeq Size:	2086 bp
RefSeq ORF:	1836 bp
Locus ID:	348807
UniProt ID:	<u>Q494V2</u>
Cytogenetics:	3q21.3
MW:	71.1 kDa
Gene Summary:	May play a role in ciliary/flagellar motility by regulating the assembly and the activity of axonemal inner dynein arm.[UniProtKB/Swiss-Prot Function]