

Product datasheet for **SC307541**

GEFT (ARHGEF25) (NM_182947) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GEFT (ARHGEF25) (NM_182947) Human Untagged Clone
Tag:	Tag Free
Symbol:	GEFT
Synonyms:	GEFT; p63RhoGEF
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC307541 representing NM_182947.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCGGGGGGGGCACAAAGGGGGTGCCTGTGCCTGTCCCCGTGTGATCCGAAAAGTGCTGGCAAAATGC
GGCTGCTGCTTCGCCCCGGGGGACGTGAATCCTATTCCATTGCGGGCAGTGAGGGGAGTATATCGGCT
TCTGCTGCCTCCGGTCTGGCTGCCCCCTTGCCCCCAGCTCTGGCCTCAGCTCTGGCCCCCTGTTCCCA
GGCCCCCAGGGCCCGTCAGTGGCCTGAGGAGATGGTTGGATCATTCAAACATTGTCTCAGTGTGGAA
ACTGAGGCAGACAGTGGTCAGGCAGGACCATATGAGAACTGGATGTTGGAGCCAGCTCTAGCCACAGGA
GAGGAGCTGCCGAACTGACCTTGCTGACCACACTGTTGGAGGGCCCTGGAGATAAGACGAGCCACCT
GAAGAGGAGACTTTGTCCAAGCCCTGAGAGTGAGGAGGAACAGAAGAAGAGCTCTGAAAGGAGT
ATGTATGCTCTGAGTGAAGTGTAGAAACAGAGAAAATGTACGTGGACGACTTGGGGCAGATTGTGGAG
GGTTATATGGCCACCATGGCTGCTCAGGGGGTCCCCGAGAGTCTTCGAGGCCGTGACAGGATTGTGTTT
GGGAATATCCAGCAAATCTATGAGTGGCAGCAGACTATTCTTGCAAGAGCTACAACGGTGTCTGAAA
GATCCTGATTGGCTGGCTCAGCTATTCAAAACACGAGCGCCGGCTGCATATGTATGTGGTGTACTGT
CAGAATAAGCCCAAGTCAGAGCATGTGGTGTGAGATTTGGGGACAGCTACTTTGAGGAGCTCCGGCAG
CAGCTGGGGCACCCTGTCAGCTGAACGACCTCCTCATAAACCTGTGCAGCGGATCATGAAATACCAG
CTGCTGCTCAAGGATTTTCTCAAGTATTACAATAGAGCTGGGATGGATACTGCAGACCTAGAGCAAGCT
GTGGAGGTCATGTGCTTTGTGCCAAGCGCTGCAACGATATGATGACGCTGGGGAGATTGCGGGGATTT
GAGGGCAAAGTACTGCTCAGGGGAAGCTCTTGGGCCAGGACACTTTCTGGGTCAACGAGCCTGAGGCT
GGAGGGCTGCTGTCTTCCCGAGGTCGAGAGAGGCGCGTCTTCTCTTTGAGCAAAATCATCATCTTCAGT
GAAGCCCTGGGAGGAGGAGTGAGAGGTGGAACACAGCCTGGATATGTATACAAGAACAGCATTAAAGTG
AGCTGCCCTGGGACTGGAGGGGAACCTCCAAGGTGACCTTGGCGCTTTGCACTGAGCTCCAGAGGGCCA
GAGGGTGGGATCCAGCGCTATGCTCTGAGGCTGCAGACCCTGCTATCAGTCAGGCTGGATCAAGCAT
GTGGCTCAGATCTTGGAGAGCCAACGGGACTTCTCAACGCATTGCAGTCACCCATTGAGTACCAGAGA
CGGGAGAGCCAGACCAACAGCCTGGGGCGGCCAAGAGGGCCTGGAGTGGGGAGCCCTGGAAGAATTCGG
CTTGGAGATCAGGCCCAGGGCAGCACACACACCCATCAATGGCTCTCTCCCTCTCTGCTGCTGTCA
CCCAAAGGGGAGGTGGCCAGAGCCCTTGGCACTGGATAAACAGGCCCTTGGTACATCCCCCAGGCT
CCCCATGACTCTCTCCAGTCTCTCAACTCAAAAACCCCTCCCTGCCAAGCCAGACTTGCCAAGCTG
GATGAAGATGAGCTGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_182947

Insert Size: 1743 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_182947.3](#)

RefSeq Size: 2563 bp

RefSeq ORF: 1743 bp

Locus ID: 115557

UniProt ID: [Q86VW2](#)

Cytogenetics: 12q13.3

MW: 63.9 kDa

Gene Summary: Rho GTPases alternate between an inactive GDP-bound state and an active GTP-bound state, and GEFs facilitate GDP/GTP exchange. This gene encodes a guanine nucleotide exchange factor (GEF) which interacts with Rho GTPases involved in contraction of vascular smooth muscles, regulation of responses to angiotensin II and lens cell differentiation. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2012]
 Transcript Variant: This variant (1) represents the longest transcript (1) and encodes isoform 1.