

## Product datasheet for SC307560

### IKZF3 (NM\_183231) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	IKZF3 (NM_183231) Human Untagged Clone
Tag:	Tag Free
Symbol:	IKZF3
Synonyms:	AIO; AIOLOS; ZNFN1A3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC307560 representing NM_183231. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG  
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC  
ATGGAAGATATACAAACAAATGCGGAAGTGAAGCACTCAGGAGCAGTCTGTGCCCGCAGAAAGTGCA  
GCGGTTTTGAATGACTACAGTTTAACCAATCTCATGAAATGGAAATGTGGACAGTGGAGAAGGCCCA  
GCCAATGAAGATGAAGACATAGGAGATGATTCAATGAAAGTGAAGATGAATACAGTGAAGAGATGAG  
AATGTTTTAAAGTCAGAACCCATGGGAAATGCAGAAGAGCCTGAAATCCCTTACAGCTATTCAAGAGAA  
TATAATGAATATGAAAACATTAAGTTGGAGAGACATGTTGTCTCATTGATAGTAGCAGGCCAACCACT  
GGAAAGATGAACTGCGATGTGTGTGGATTATCCTGCATCAGCTTCAATGTCTTAATGGTTCATAAGCGA  
AGCCATACTGCAAGTGCGGAGGCAAGACACATCAAAGCAGAGATGGGAAGTGAAGAGCTCTCGTACTG  
GACAGATTAGCAAGCAATGTGGCAAAACGAAAAAGCTCAATGCCTCAGAAATTCATTGGTGAGAAGCGC  
CACTGCTTTGATGTCAACTATAATTCAAGTTACATGTATGAGAAAGAGAGTGAGCTCATACAGACCCGC  
ATGATGGACCAAGCCATCAATAACGCCATCAGCTATCTTGGCGCCGAAGCCCTGCGCCCTTGGTCCAG  
ACACCGCTGCTCCACCTCGGAGATGGTTCAGTTATCAGCAGCATGTATCCATAGCCCTCACCCGG  
GCTGAGATGTCAAACGGTGCCCTCAAGAGCTGGAAAAGAAAAGCATCCACCTCCAGAGAAGAGCGTG  
CCTTCTGAGAGAGGCCTCTCTCCAACAATAGTGGCCACGACTCCACGGACACTGACAGCAACCATGAA  
GAACGCCAGAATCACATCTATCAGCAAAATCACATGGTCTGTCTCGGGCCGCAATGGGATGCCACTT  
CTGAAGGAGGTTCCCGCTCTTACGAACCTCTCAAGCCCCGCCCCATCTGCCAAGAGACTCCGTCAAA  
GTGATCAACAAGGAAGGGGAGGTGATGGATGTGTATCGGTGTGACCACTGCCGCGCTCTCTCTGGAC  
TATGTGATGTTACGATTACATGGGCTGCCACGGCTTCCGTGACCCCTTCGAGTGAACATGTGTGGA  
TATCGAAGCCATGATCGGTATGAGTTCTCGTCTCACATAGCCAGAGGAGAACACAGAGCCCTGCTGAAG  
TGA  
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT  
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC



View online »

<b>Restriction Sites:</b>	Sgfl-Mlul
<b>ACCN:</b>	NM_183231
<b>Insert Size:</b>	1245 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_183231.2</a></u>
<b>RefSeq Size:</b>	9401 bp
<b>RefSeq ORF:</b>	1245 bp
<b>Locus ID:</b>	22806
<b>UniProt ID:</b>	<u><a href="#">Q9UKT9</a></u>
<b>Cytogenetics:</b>	17q12-q21.1
<b>MW:</b>	47 kDa

**Gene Summary:**

This gene encodes a member of the Ikaros family of zinc-finger proteins. Three members of this protein family (Ikaros, Aiolos and Helios) are hematopoietic-specific transcription factors involved in the regulation of lymphocyte development. This gene product is a transcription factor that is important in the regulation of B lymphocyte proliferation and differentiation. Both Ikaros and Aiolos can participate in chromatin remodeling. Regulation of gene expression in B lymphocytes by Aiolos is complex as it appears to require the sequential formation of Ikaros homodimers, Ikaros/Aiolos heterodimers, and Aiolos homodimers. Several alternative transcripts encoding different isoforms have been described, as well as some non-protein coding variants. [provided by RefSeq, Apr 2012]

Transcript Variant: This variant (5) lacks two alternate in-frame exons compared to variant 1. The resulting isoform (5, also known as Aio-del4,5) has the same N- and C-termini but is shorter compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.