

Product datasheet for **SC307706**

CFAP20DC (NM_198463) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CFAP20DC (NM_198463) Human Untagged Clone
Tag:	Tag Free
Symbol:	CFAP20DC
Synonyms:	C3orf67
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC307706 representing NM_198463.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGATCAAACGTAAAATTTGGTGCATCTATGCATTGACTTAGTAGCATTACCAGTGAATATTTCAAG
GGGGCAGTTTTCCAGTCATTGGATGGAATTGTTGTCTCAGCTAACTGTAAGCTACGGAAGATCTTCACC
TTAAAATCAAAGCCACAAGACACTGCTGATAAGGATGCTGTCTATGGTGTCCCTTTTCAACAGATGAG
CCTACAGATATTATACCAGGAAGCTGTCAACTAATGACAGATGTTCCACATGTCACACAGCTGCTAAAC
ATGACTAAACTTCGCCAAACTGAAATAAAATTCGGAGGCCATCCTCTAAGATCAGCAGAATCAGATCAG
TTCATTAACAGAGGAACAAGTATTACACGGAACAGTAAAATCAAGATGTTTGTCAATCGCATTGGGA
TCCAAAGTTCTTGGGCCACCTCCACTCTCTGGCAGAAGGAATAACATGAAGATATCCAGCGAGACAGTG
AGATCCGTTGGGTCAAAAATAACCGATCATGCCAGCCGCTCCACTGTAGAGAAGTGTGTTAATGGTACA
GAAATGTCAGCCTTGCTGATACCTGAGTCTGAGGAACAAGGAAAATAAAGAAAATATTCACCAATAAAG
CAGACTGTACCTATTCATGCAGCCAATCTACATATTATGCATCCGCATCCCCCTCAAGAACCATCAGCA
GATAAGAATAATAACAGAAGAAGATTACGGTAAAAAGTACCAGCAGAGAAGGACAGAGACACCCAGC
GGTAGCTCTTCAGGAAATAATAGGATTGAAGATAAAGCATCAACTATCCTCACCCTGTGTCCCAACAA
GGAGCAGAGCTGTTGAACTCCGGCACTCTAGGACCCAGTCTCCTGATCAATCAGATGAGTGGATTTTT
CCTGAAAATGCTGATCACATTTTCATATCTGGCATCCAGCAGACAGTCTCTACTTCTGGGTGATGACTCC
TGCAACCCATCACACCTGTGGCTGGAAGCCAGCAAAGAGAGTGAACACGACCAGCAGGCAGAGGAATCC
CAGAGTGTCCAAAGGACATTTTCACTTTTTCATCAAGACCACGATCAGCACCTCATGAAAAGACTCAG
ACTATGTCCCCAGAGGAGCTCTCATTTATTTGGATCTAAAAGAGGATAACAGTGTGACAAGCAGAGAC
ACCCAATCAGAGGATGATTTTTACGGCGGCAGCAGCAGTGAAGAGGAATGACTGGCGAAACTATCAG
CCAAGCCAGATGAGTGAATCCGAGTTACAGATGCTAGCAAGCCTACGGTGGCAACAAAATGAAGAATG
GAGGATGCTGGGACCTCCCATGGCTGAGTGCCTCCCAGGTGGACAACCTGTAATGTCAGCATAAGTACC
AGCAGTGACGACACAACCACCTGGAACCTCTGCCTGCCACCCCTGTCAACCAGGGTCGCCACTATCAG
AAAGAAATGAACCCACCTTCTCCTTCTAATCCCGGGACTGGTTAAATATGTTGAGCCACCAATCGTT
CCTCCCAGTCAACAGCCGGCTGAGCAGCGTCCAGATTCTGTGAAAGTTTGAGTGTCAAGGTGAAGAA
GACCTCAGTGTGGAAGAGGACGAGGAAGTACTGACTTTGTTGTATGACCCTTGTCTGAACGTACTTT
GACCCCAAACAGGAAATACTATGAGTTGGTATAA
ACGGGTACGGCGCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_198463

Insert Size: 1692 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_198463.2](#)

RefSeq Size: 2551 bp

RefSeq ORF: 1692 bp

Locus ID: 200844

UniProt ID: [Q6ZVT6](#)

Cytogenetics: 3p14.2

MW: 62.8 kDa