

Product datasheet for SC307715

VSTM1 (NM_198481) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	VSTM1 (NM_198481) Human Untagged Clone
Tag:	Tag Free
Symbol:	VSTM1
Synonyms:	SIRL-1; SIRL1; UNQ3033
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC307715 representing NM_198481. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACCGCAGAATTCCTCTCCCTGCTTTGCCTCGGGCTGTGTCTGGGCTACGAAGATGAGAAAAAGAAT
GAGAAACCGCCCAAGCCCTCCCTCCACGCCTGGCCAGCTCGGTGGTTGAAGCCGAGAGCAATGTGACC
CTGAAGTGTGAGGCTCATTCCAGAAATGTGACATTTGTGCTGCGCAAGGTGAACGACTCTGGGTACAAG
CAGGAACAGAGCTCGGCAGAAAACGAAGCTGAATCCCTTCACGGACCTGAAGCCTAAGGATGCTGGG
AGGTACTTTTGTGCTACAAGACAACAGCCTCCCATGAGTGGTCAGAAAGCAGTGAACACTTGCAGCTG
GTGGTCACAGATAAACACGATGAACTTGAAGCTCCCTCAATGAAAACAGACACCAGAACCATCTTTGTC
GCCATCTTCAGCTGCATCTCCATCCTTCTCCTCTCCTCTCAGTCTTCATCATCTACAGATGCAGCCAG
CACAGTTCATCATCTGAGGAATCCACCAAGAGAACCAGCCATTCCAAACTCCGGAGCAGGAGGCTGCC
GAGGCAGATTTATCCAATATGAAAGGGTATCTCTCTCGACGGCAGACCCCAAGGAGTGACCTATGCT
GAGCTAAGCACCAGCGCCCTGTCTGAGGCGCTTCAGACACCACCAGGAGCCCCAGGATCTCATGAA
TATGCGGCACTGAAAGTGTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_198481
Insert Size:	711 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_198481.3</u>
RefSeq Size:	1034 bp
RefSeq ORF:	711 bp
Locus ID:	284415
UniProt ID:	<u>Q6UX27</u>
Cytogenetics:	19q13.42
Protein Families:	Transmembrane
MW:	26.1 kDa
Gene Summary:	<p>Isoform 2: Behaves as a cytokine, promoting IL17A secretion by CD4+ T-cells, and differentiation and activation of IL17 producing helper T-cells (TH17).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) encodes the longest isoform (1).</p>