

Product datasheet for SC307717

ZNF621 (NM_198484) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF621 (NM_198484) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF621
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC307717 representing NM_198484. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTCCAAACAACCTGGCCTCAGGAGTCAGTGACCTTTGAGGATGTGGCTGTTACTTCACCCAGAAT
CAATGGGCCAGCCTCGACCCTGCGCAGAGGGCCCTGTACGGGGAGGTGATGCTGGAGAATTATGCAAAT
GTGGCTTCTCTGGTAGCGTTTCCATTCCCCAAACCTGCTGATCTCCCACCTGGAGAGAGGGGAAGCA
CCATGGGGCCAGATCCCTGGGACCCGAGATTCTGAGAGGCATCAGTCAAGGTGGTGGTCTGGATC
AAAAATGAAGGGCTAGTTATAAAGCAGGAAGCCTCTGAAGAAACAGAGTTGCACAGAATGCCAGTAGGA
GGACTTCTCAGGAACGTTTCTCAGCACTTTGATTTTAAAGGAAGGCCTGAAGCAGACTTCAATCTA
AATCCAAATCTGATACTTCGAGGTGGAATGAAGTTCTATGAATGTAAAGAATGTGGGAAAATCTCCGA
TATAACTCAAAGCTTATTCGGCATCAGATGAGTCATACTGGGAAAAGCCCTTAAAGTGAAGGAGTGT
GGCAAAGCTTCAAGTCCAGCTATGATTGTATTGTACATGAGAAAAACCACATTGGAGAAGGGCCCTAT
GAATGTAAGGAGTGTGGCAAAGTTTGGATTCCAACACAGCCTTGACTCAACATCAGAGGATCCACACT
GGAGAGAAACCCTATGAATGTAAGAGTGTGGAAAGGCTTCCGTAGGAGTGCAGCATACCTGCAGCAT
CAGAGATTACACACGGGAGAGAACTCTATAAATGTAAGGAATGTTGGAAAGCTTTCGGTGTAGGTCA
CTTTTATTGTCCATCAGAGAATTCATACTGGGAGAAACCTTATCAATGTAAGGAGTGTGGCAAAGCC
TTCACCCAGAAAATAGCCTCCATTCAGCATCAGAGAGTTCACACTGGGGAGAAGCCTTATGAATGTAAG
GTGTGTGGGAAAGCCTTCAAATGGTATGGAAGTTTGTTCAGCATCAGAAATGCACCCTGTGGAGAAG
AAGCCAGTCAAGGTCCTTGGGCCATCCCTGGTCACTCCCAAGTGTCTCTCCAGCCATACCTCCTGTT
CTTCTCCAGGGATCCTGTTCTGCTTCAGCCGTAGCTGTGCCTTCACTGACCTTCCACATGCTGTGCTC
ATTCTACCTCTGGGAATTTTTTCATGCTGCTGCCTACATCTGGAATACCTTCTTCATCTGCCAAATA
GTGCGTGTCTCCAGGTCTTACTCCACTGTGAAACCTCCCAAGTATTCTCACCCCTTCTTCTCAC
TCCTCATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

ACCN:	NM_198484
Insert Size:	1320 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_198484.4
RefSeq Size:	8483 bp
RefSeq ORF:	1320 bp
Locus ID:	285268
UniProt ID:	Q6ZSS3
Cytogenetics:	3p22.1
MW:	49.2 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longest transcript and encodes the longer isoform (1). Both variants 1 and 2 encode the same isoform. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.