

Product datasheet for SC307718

TPRG1 (NM_198485) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TPRG1 (NM_198485) Human Untagged Clone
Tag:	Tag Free
Symbol:	TPRG1
Synonyms:	FAM79B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC307718 representing NM_198485. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
 ATGTCAACAATTGGGAGTTTGAAGATTCCAGGCTGTGTCTCTGAAGCAAGAGGGAGATGACCAACCC
 TCTGAGACTGACCACCTATCGATGGAGGAAGAGGACCGATGCCAAGACAGATTTCAGGCAGTCAAGT
 GTGACCGAATCAACTCTTACCCCAATCCTTATCATCAGCCTTATATCTCACGGAAGTACTTTGCTACA
 CGGCCGGGGGCCATTGAGACTGCCATGGAAGACTTGAAAGGTACGTAGCTGAGACTTCTGGAGAGACC
 ATCAAGGCTTCTGGCTCTTGACAAAGATAGACCACTGGAACAATGAGAAGGAGAGAATTCTACTGGTC
 ACAGACAAGACTCTTGTATCTGCAAATACGACTTCATCATGCTGAGTTGTGTGCAGCTGCAGCGGATT
 CCTCTGAGCGCTGTCTATCGCATCTGCCTGGGCAAGTTCACCTTCCTGGGATGTCCCTGGACAAGAGA
 CAAGGAGAAGGCCTTAGGATCTACTGGGGAGTCCGGAGGAGCAGTCTCTTCTGTCCCGCTGGAACCCA
 TGGTCCACTGAAGTTCCTTATGCTACTTCACTGAGCATCCTATGAAATACACCAAGTGAAGAAATCCTT
 GAAATTTGCAAGTTGTCTGGGTTTCAAGCTTGTTCAGCTATCCAGAATGCCCACAAGAATTCA
 ACTGGATCTGGAAGAGGAAAGAACTGATGGTGTAACTGAACCCATTTTGATTGAGACCTACACAGGG
 CTGATGTCATTCATTGGAACCGCAACAACTTGCTATCCCTTGCCCGTGGGAGTATTGGTTTTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_198485
Insert Size:	828 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_198485.3</u>
RefSeq Size:	3665 bp
RefSeq ORF:	828 bp
Locus ID:	285386
UniProt ID:	<u>Q6ZUI0</u>
Cytogenetics:	3q28
MW:	31.2 kDa