

Product datasheet for **SC307743**

B3GNT8 (NM_198540) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	B3GNT8 (NM_198540) Human Untagged Clone
Tag:	Tag Free
Symbol:	B3GNT8
Synonyms:	B3GALT7; beta3Gn-T8; BGALT15
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC307743 representing NM_198540. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCGCTGCCCAAGTGCCTTCTCTGCCTGTCAGCACTGCTCACACTCCTGGGCCTCAAAGTGTACATC
GAGTGGACATCCGAGTCCCGGCTCAGCAAGGCCTACCCAGCCCTCGGGGCACCCCGCAAGCCCCACG
CCAGCCAAACCTGAGCCACCCCTACCTGCCAACCTCTCCACCCGCTGGGCCAGACTATCCCGCTGCC
TTTGCTTACTGGAACCAGCAGCAGTGGCGGCTGGGGTCCCTGCCAGTGGGGACAGCACTGAAACGGGG
GGCTGCCAGGCTTGGGGGGCCCGCCGCCACCGAGATCCCTGACTTCGCCTCCTACCCCAAGGACCTC
CGCGCTTCTTGCTGTCAGCAGCCTGCCGAGCTTCCACAGTGGCTGCCTGGAGGTGGTGGCAGCCAA
GTCTCCAGCTGCTCAGATACTGATGTCCCTACCTGCTGTTGGCCGTCAGTCAAGTCAAGACAGGGCGCTT
GCAGAACGACAGGCCGTGAGAGAGAGCTGGGGCAGTCCAGCTCCAGGGATCCGGCTGCTCTTCTGCTA
GGGTCTCCGGTGGGTGAGGCGGGGCTGACCTAGACTCACTAGTGGCTGGGAGAGCCGTCGCTACAGT
GACCTGTGCTCTGGGACTTCTCGACGTCCCATCAACCAGACGCTCAAAGACCTGCTGCTGCTGGCC
TGGCTGGGCCGCCACTGCCCCACCGTAGTTTTGTCTTGCGAGCTCAGGACGATGCCTTTGTACACACC
CCTGCCCTGCTGGCTCACCTGCCGGCCCTGCCACCTGCCTCGGCCGAAGCCTCTACCTGGGTGAGGTG
TTTACCAGGCCATGCCTCTCCGAAGCCAGGAGGACCTTCTATGTGCCGAGTCTTCTCGAAGGT
GGCTACCAGCCTATGCAAGCGGGGTGGCTACGTCATTGCCGGGCGCTGGCACCTGGCTGCTGCGG
GGCGCAGCCCGTGTGGCACCTTCCCCTTTGAGGACGCTCACTGGCCCTTGCATCCGAGCCCTGGGC
CTGGTGCCCCAGGCCACCCAGGCTTCTCACAGCCTGGCCAGCAGACCCGACTGCGGACCACTGTGCT
TTCCGCAACCTGCTGCTGGTACGGCCCTGGGCCCCAGGCCAGCATTGGCTCTGGAAACAACCTGCAA
GACCAAGGCTCCAGTGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul



[View online >](#)

ACCN:	NM_198540
Insert Size:	1194 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_198540.2
RefSeq Size:	1894 bp
RefSeq ORF:	1194 bp
Locus ID:	374907
UniProt ID:	Q7Z7M8
Cytogenetics:	19q13.2
MW:	43.4 kDa
Gene Summary:	Beta-1,3-N-acetylglucosaminyltransferase that plays a role in the elongation of specific branch structures of multiantennary N-glycans. Has strong activity towards tetraantennary N-glycans and 2,6 triantennary glycans.[UniProtKB/Swiss-Prot Function]