

Product datasheet for **SC307746**

SPATA21 (NM_198546) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SPATA21 (NM_198546) Human Untagged Clone
Tag:	Tag Free
Symbol:	SPATA21
Synonyms:	spergen-2; spergen2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC307746 representing NM_198546.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGACAATAGAAACACCCAGATGTACACGGAGGAGGAAAAGACAGTCAACCCCTTCCTGCCATCCAGC
CCTGGACCTAAAAAGCAAAGGAGGGGTGAGGCTGTGGAGACCCACCCAGCACCGGGGCTCTTCCT
CCACCAGAGGTGAGGGACATTGGAGAGAGGCGGGAGCCAGACCGTGCAGCAGCAGCCTCAGAAGCCC
GCGGTGGCTGCAGGGACACAGAGCCTCGGGAACCTCCGGCAGGGCTTCATGAAGTCTTGTGGAGGTG
GAGAAAATGGAGGCCTCCCATCGGAGAGCCTCGAAGGCCCGGTCCAGACAGCCAGAAGTCGCCAGG
ACCCTGACCCCGGTGCCACCTCAGCTCCAAGCCTGCCTCAGACCCCTGCCTCGGTCCCTGCCAGCGGC
CCATCGTGGGCCCGGTGCCAGCTCCTGGGCCAGAACCTGCCCCATGGGAGCCCGGTCCACCTCC
ATGCTTGGCCTGTCTGCTGGGCCCTGCCCTGGACCTGGGTGGAGAAGGATGGAATCTTGACACAG
AGCAGCGAGAGAACCCTGAGCTACGCCAAGGCGCGGCAGGAGCCGAAGAGCAGAGCCTCCAAAAGCTT
TATCAAAACCGGAGAGAAGTCCGAGGAGCAACTGACCCTGAAGCAAGAGGAAGCCTCCGCAGCTACTTT
GAGATCTTCAATGGTCTCGGTGAGGTGGATGCACAGAGCCTGAAGAATATCCTGCTCCTAATGGGCTTC
TCTGTGACGCTGGCCAGGTGGAGGACGCCCTGATGAGTCTGATGTCAATGGAGATGGTCTGTGGAC
TTCAAAGACTTCTTGGCTGTGATGACAGACACAGGCGCTTCTTCTGCTCTGTGGAACAGAACGCCCTG
TCGGACATGGCTCCCCACAACCCCCACACTCTACTCTTTGAGATCCTGTCCCTGCTGGTAGAGATGCTG
GCCTTACCAGAGGCAGTCTTAGAGGAAATCACAACTACTACCAGAAAAAGCTGAAGGAAGGCACCTGC
AAGGCCAGGAGATGGAAGCGGCCGTAGGCCGGCTGCGGTTGCAGAAGCTTCCCTACAACCCCCAGCAA
GAAGAGAGCTCAGAAGTCCAGAACGAAAGTCTCAGTATCCTGAGCCGGCTGAAGCAGCAGAACTAT
GCTCCCAACCTGCAGAGCCCTATGCCAGGTGCCCTGCATCCTGCTCTGCCACAGCTGGACAAGAAG
ATGGTCCGTAGGCAGCCGAGCAACCACTATGCACTAGACCAGTGCACACCCCTGGCCTGGATCCTGAC
ATCCGCAGCCCTTCTTCCAGTCAGGATCTCAAGGAAACAGGGAACACAACCTCTGACAGCAGAAAGTGG
CTCAGCTCTGTGCCAGCTCGAACCCACTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_198546

Insert Size: 1410 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_198546.1</u>
RefSeq Size:	2015 bp
RefSeq ORF:	1410 bp
Locus ID:	374955
UniProt ID:	<u>Q7Z572</u>
Cytogenetics:	1p36.13
MW:	52.2 kDa
Gene Summary:	Involved in the differentiation of haploid spermatids.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longer isoform (1).