

Product datasheet for **SC307766**

KRTAP10-4 (NM_198687) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KRTAP10-4 (NM_198687) Human Untagged Clone
Tag:	Tag Free
Symbol:	KRTAP10-4
Synonyms:	KAP10.4; KAP18-4; KRTAP10.4; KRTAP18-4; KRTAP18.4
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_198687, the custom clone sequence may differ by one or more nucleotides <pre> ATGTCCGTCTGCTCCAGCGACCTGAGCTACAGCAGCCGCTGCTGCCTTCCTGGTTCTGT GACTCTTGCTCCGACTCCTGGCAGGTGGACGACTGCCAGAGAGCTGCTGCGAGCCCCC TGCTGCGCCCCCAGCTGCTGCGCCCCGCCCCCTGCCTGAGCCTGGTCTGCACCCAGTG AGCCGTGTGTCCAGCCCCCTGCTGCCAGTGACCTGCGAGCCCAGCCCCTGCCAATCAGGC TGACACGACTCCTGCACGCCCTCGTGCTGCCAGCAGTCTAGCTGCCAGCTGGCTTGCTGT GCCTCCTCCCCCTGCCAGCAGGCCTGCTGCGTGCCCGTCTGCTGCAAGACTGTCTGCTGC AAGCCTGTGTGTGTGTGCCCCGTCTGCTGTGGGGATTCTTCATGTGCCAGCAGTCTAGC TGCCAGTCAGCTTGCTGCACCTCCTCCCCCTGCCAGCAGGCCTGCTGTGTGCCCATCTGC TGCAAGCCTGTCTGCTCTGGGATTTCCTCTTCGTGCTGCCAGCAGTCTAGCTGTGTGAGC TGTGTGTCCAGCCCCTGCTGCCAGGCGGTCTGTGAGCCCAGCCCCTGCCAATCAGGCTGC ATCAGCTCCTGCACGCCCTCGTGCTGCCAGCAGTCTAGCTGCCAGCCGGCTTGCTGCACC TCCTCCTCCTGCCAGCAGGCCTGCTGCGTGCCCGTCTGCTGCAAGACTGTCTGCTGCAAG CCTGTGTGCTCTGAGGATTCTCTTCATGTGCCAGCAGTCTAGCTGCCAGCCGGCTTGCTG TGACCTCCTCTCCCTGCCAGCAGGCTTGCTGTGTGCCTGTCTGCTGCAAGCCTGTGTGC TGCAAGCCTGTGCGCTCTGTGCCATCTGCTCTGGGGCTTCCTCTCTGTGCTGCCAGCAG TCTAGCTGCCAGCCAGCTTGCTGCACCTCCTCCCAAAGCCAGCAGGGCTGCTGCGTGCCC GTCTGCTGCAAGCCTGTGAGCTGTGTGCCTGTTTGTCTGCGGGCTTCCTCTTCATGCTGC CAGCAATCTAGCTGCCAGCCAGCTTGCTGCACCACTCCTGCTGCAGACCCTCCTCCTCC GTGTCCCTCCTTGCCGCCCCGTGTGCAGGCCCGCTGCTGCGTGCCCGTCCCTTCCTGC TGTGCTCCACCTCCTCTGCCAACCCAGCTGCTGCCGCCAGCCTCCTGCGTGTCCCTC CTCTGA </pre>
Restriction Sites:	Please inquire
ACCN:	NM_198687
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_198687.1</u> , <u>NP_941960.1</u>
RefSeq Size:	1382 bp
RefSeq ORF:	1206 bp
Locus ID:	386672
Cytogenetics:	21q22.3
Gene Summary:	This is an intronless gene located in a cluster of related genes on the q arm of chromosome 21. The proteins encoded by these genes form disulfide bonds with cysteine residues in hair keratins, thereby contributing to the structure and stability of hair fibers. [provided by RefSeq, Apr 2014]