

Product datasheet for SC307820

OR7C1 (NM_198944) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: OR7C1 (NM_198944) Human Untagged Clone
Tag: Tag Free
Symbol: OR7C1
Synonyms: CIT-HSP-146E8; HSTPCR86P; OR7C4; OR19-5; TPCR86
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >SC307820 representing NM_198944.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAACAGGAAATCAAACACATGCCAAGAATTTCTCCTCCTGGGATTTTCAGCAACGTCAGAGATT
CAGTTCATTCTCTTTGGGCTGTTCTCTCCATGTACCTAGTCACTTTACCAGGAACTGCTCATCATC
CTGGCCATATGCTCAGACTCCCACCTCCACACCCCATGTACTTCTTCTCTCCAACCTGTCTTTTGCT
GACCTCTGTTTTACCTCCAGACTGTCCAAAGATGTTACTGAATATACTGACACAGAACAAATTCATA
ACATATGCAGGCTGTCTCAGTCAGATTTTTTTTTTCACTTCATTTGGATGCCTGGACAATTTACTCTTG
ACCGTGATGGCCTATGACCGCTTCGTGGCGTCTGTCAACCCCTGCACTATACGGTCATCATGAACCCC
CAGCTCTGTGGACTGCTGGTTCTGGGGTCTGGTGCATCAGTGCATGGGTTCCCTGCTCGAGACCTTG
ACTGTTTTGAGGCTGTCTTCTGCACCGAAATGGAAATCCACACTTTTTTTGTGATCTACTTGAAGTC
CTGAAGCTCGCCTGTTCTGACACCTTCAATAAACGTGGTGATATACTTTGCAACTGGCGTCTGGGT
GTGATTTCTTCACTGGAATATTTTTCTTACTATAAAATGTTTTCTATACTGAGGATTTCTCA
GCTGGGAGAAAGCACAAGCGTTTTCCACCTGTGGTCCCACCTCTCAGTGGTCACTTGTCTATGGC
ACGGGCTTTGGGGTCTATCTCAGTCTGCAGCCACACCATCTTAGGACAAGTCTGGTGGCCTCAGTG
ATGTACACCATGGTCACCCCATGTGAACCCCTTCATCTACAGCCTGAGGAACACGGACATGAAGAGG
GCCCTGGGAGACTCCTCAGTAGGGCAACATTTTTTAATGGTGACATCACTGCAGGACTTTCAATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI
ACCN: NM_198944
Insert Size: 963 bp



[View online »](#)

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_198944.1](#)

RefSeq Size: 963 bp

RefSeq ORF: 963 bp

Locus ID: 26664

UniProt ID: [O76099](#)

Cytogenetics: 19p13.12

Protein Families: Transmembrane

Protein Pathways: Olfactory transduction

MW: 35.5 kDa

Gene Summary:

Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]