

Product datasheet for **SC307869**

ZNF585A (NM_199126) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF585A (NM_199126) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF585A
Synonyms:	Zfp27
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF:	>NCBI ORF sequence for NM_199126, the custom clone sequence may differ by one or more nucleotides <pre> ATGCTGGAGACCTACAGCCACCTGCTCTCAGTAGGATATCAAGTTCCTGAAGCAGAGGTG GTCATGTTGGAGCAAGGAAAGGAACCATGGGCACTGCAGGGTGAGAGGCCACGTGAGAGC TGCCCAGGAGAGAAATATGGGACCATAATCAATGTAGAAAAATCCTCAGTTATAAACAA GTATCCTCTCAACCTCAAAAAATGTATCCTGGGAGAAAGCTTATGAATGCGCCAAATTT GAAAAGATATTCACCCAGAAGTCACAGCTCAAAGTACACCTGAAAGTTCTTGAGGAGAA AAGCTCTATGTATGCATTGAATGTGGGAAGGCTTTGTACAGAAGCCAGAATTTATTATA CACCAGAAAACCCATATGAGAGAGAAACCTTTAAATGCAATGAATGTGGAAAATCCTTT TTTCAAGTGTCGTCCTCTTCAGGCATCAGAGAATTCATACCGGAGAGAACTCTATGAA TGCAGCCAGTGTGGGAAGGCTTCTCTTATAACTCAGATCTCAGTATACATGAGAAAATT CATACTGGAGAGAGACCATGAATGCACTGACTGTGGCAAAGCATTACACAAAAGTCC ACACTCAAGATGCATCAGAAAATCCATACAGGCGAGAGATCTACATCTGTATTGAATGC GGACAGGCCTTCATCCAGAAGACCCATTTGATTGCACACCGAAGAATTCATACTGGAGAA AAACCATATGAGTGCAGTAAGTGTGGCAAATCCTTCATTTCCAAGTCACAACTTCAGGTA CATCAACGTGTTACACAAAGAGTGAAGCCCTATATATGTACCGAATATGGGAAGGCTTTC AGCAATAATTCCAACCTCGTTACACATAAGAAAGTTCAAAGTAGAGAGAAATCTTCCATA TGTAAGTGTGGGAAGGCTTTACCTACAGGTCAGAGTTGATTATTCATCAGAGAATT CACACTGGAGAGAAACCTTATGAATGCACTGACTGTGGGAAGGCTTCACTCAGAAGTCA GCACTCAGAGTGCATCAGAGAATTCATACAGGAGAAAAATCGTATATATGCATGAATGT GGACTGGCCTTCATTCAGAAGGCACACTTGATTGCACATCAAATAATTCATACTGGAGAG AAACCTCATAAATGTGGTCACTGTGGGAATTTGTTTACCTCCAAGTCGCAACTCCATGTT CATAAACGAATTCACACAGGAGAAAAAGCCCTATATGTGCAATAAATGTGGGAAGGCATTTC ACCAACCGGTCAAATCTCATTACACATCAGAAAACCTCATACAGGAGAGAAATCTTATATA TGTTCCAATGTGGGAAGGCTTACCCAGAGGTGAGACTTGATTACACATCAGAGAATC CATACTGGGAGAGGCTTATGAATGCAATACTTGTTGGAAAAGCCTTCACTCAGAAGTCA CACCTCAATATACATCAGAAAATTCACACTGGAGAGAGACAGTATGAATGCCACGAATGT GGGAAAGCCTTCAACCAGAAATCAATACTCATTGTTTCATCAGAAAATTCATACAGGAGAG AAACCTATGTATGCACTGAGTGTGGAAGAGCTTTCATCCGAAGTCAAACCTTTATTACT CATCAAAGAATTCATACTGGAGAGAAAGCCTTATGAATGCACTGACTGTGGGAAGTCTTTT ACCTCCAAGTCTCAGCTCCTGGTGCATCAGCCAGTTCACACAGGAGAGAAACCTATGTG TGTGCCGAGTGCGGAAGGCTTTAGTGGCAGGTCAAATCTCAGTAAGCACCAGAAAACT CATACCGGAGAAAAGCCCTACATTTGTTCTGAATGTGGGAAGACCTTTCGACAGAAGTCA GAGTTGATTACACATCAGAGAATTCATACTGGAGAGAAACCTTATGAGTGCAGTGAAGT GGGAAGTCTTTCATAAAAAATCACAGCTCCAAGTGCATCAGCGAATTCACACTGGAGAG AAGCCTTACGTGTGTGCTGAGTGTGGGAAGGCTTTACTGACAGGTCCAATTTGAATAAA CATCAGACAACACACTGGAGACAAACCTACAAGTGTGGCATCTGTGGGAAGGCTTTC GTTCAGAAATCAGTGTTACGCTCCATCAGAGCAGCCACGCTTGA </pre>
Restriction Sites:	Please inquire
ACCN:	NM_199126
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_199126.1, NP_954577.1</u>
RefSeq Size:	4107 bp
RefSeq ORF:	2145 bp
Locus ID:	199704
UniProt ID:	<u>Q6P3V2</u>
Cytogenetics:	19q13.12
Protein Families:	Transcription Factors
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) has an alternate 5' exon and an additional exon in the 5' region, which results in translation initiation at a downstream AUG, compared to variant 3. The resulting isoform (a) has a shorter N-terminus, compared to isoform b. Both variants 1 and 2 encode the same isoform a. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.