

Product datasheet for **SC307874**

FOXD4L2 (FOXD4L3) (NM_199135) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FOXD4L2 (FOXD4L3) (NM_199135) Human Untagged Clone
Tag:	Tag Free
Symbol:	FOXD4L3
Synonyms:	FOXD4B; FOXD4L2; FOXD4L4; FOXD6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC307874 representing NM_199135. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTAACCGTCAGAATTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAACCTGCCAAGAGCTGAGCGCCTTCGCTCCACACCGCAGCGCAGCCTCCGGGACTCCGATGGGGAA
GACGGTAAAATCGATGTCCTGGGAGAGGAGGAAGATGAAGACGAGGTGGAAGACGAGGAGGAGCGCGC
AGCCAGCAGTTCCCTAGAGCAGTCGCTCCAGCCGGGGCTGCAGGTGGCCCGGTGGGCGGGGTTGCGCTT
CCCCGAGAGCACATCGAGGGCGGGCGGGCCCGAGCGACCCCTCAGAGTTTGGCACCAAGTTCAGGGCA
CCGCCAAGGTCTCGGCGGCCTCTGAAGATGCCCGGCAGCCGGCAAAGCCCCCTACTCGTACATCGCG
CTCATCACCATGGCCATCTGCAAAACCCGCACAAGCGCCTCACGCTCAGCGGCATCTGCGCCTTATT
AGTGGCCGCTTCCCTACTACCGCCGCAAGTTCGCCGCTGGCAGAACAGCATCCGCCACAACCTCTCG
CTGAACGACTGCTTCTGTTAAGATCCCCCGCAGCCGGGCCACCCAGGCAAGGGCAACTACTGGAGCCTG
GACCCCGCTCCCAAGACATGTTTCGACAATGGCAGCTTCTCCGGCGTAGGAAGCGTTTCAAGCGCCAC
CAACTGACCCCGGAGGCCACCTGCCCCACCCCTTCCCTACTCTGCTGCACACGCGCCCTGCACAAC
CCCCGCCAGGCCCTCTGCTTGGGGCCCTGCCCGCCGAGCCAGTCCCGGGGGCTACCCCAACACC
GCCCCCGGAGACGCCCTTACGCTCTGCTGCACCCGCATCCTCTTCGCTACTACTGCTCTCGGCCCCC
GTCTATGCCGGGGCACCAAGAAAGCAGAAGGCGGGCCCTGGCGTCCGACCCGGCACCCCTCCCGTGTGC
AGCCCTCACTTGGTCTCAGCCTTGGGAGGAGGCAAGGTCTGGCGTCCGACCCGGGAGGCGGATGCA
TCTTTTCAGCATTGAGAGTATTATGCAAGGGGTGAGGGGAGCGGTACAGGGGCTGCGCAGAATTTGT
CCCCGACCGCGTGGAGCTACTGCCACCTGCTCCAGCGACCATCAAGCCTGTTGCATCCCCAGACCGCTG
CCCTTTGCTGCAAGTGTCCGCCGCCGCTGCTCGGACAATTTGCAGCAATAGCAGCAGCATCAGG
AGGAGGACTGCGCCAACGGCTGCGCTCCACCAAGGGCGCGGTGCTGGGCGGGCACCTGTGCGCCTCGT
CGGCCCTGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online >](#)

Restriction Sites:	Sgfl-Mlul
Plasmid Map:	□
ACCN:	NM_199135
Insert Size:	1254 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_199135.4</u>
RefSeq Size:	2218 bp
RefSeq ORF:	1254 bp
Locus ID:	286380
UniProt ID:	<u>Q6VB84</u>
Cytogenetics:	9q21.11
MW:	45.8 kDa