

## Product datasheet for **SC307901**

### **TCL1B (NM\_199206) Human Untagged Clone**

#### Product data:

Product Type:	Expression Plasmids
Product Name:	TCL1B (NM_199206) Human Untagged Clone
Tag:	Tag Free
Symbol:	TCL1B
Synonyms:	T-cell leukemia/lymphoma 1B; T-cell lymphoma/leukemia 1B; TCL1/ MTCP1-like 1; TML1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC307901 representing NM_199206. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG  
 GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**  
 ATGGCCTCCGAAGCTTCTGTGCGTCTAGGGGTGCCCTGGCGTCTGTGGATCCAGAGGCCTGGCATC  
 TACGAAGATGAGGAGGGGAGAACCTGGGTGACTGTGGTCGTGCGTTCAATCCCTCGCGTAGGGAATGG  
 GCCAGGGCCTCCAGGGCAGCAGATATGAACCCAGCATCACAGTGCACTTGTGGCAGATGGCAGTGCAT  
 ACCCGGGAGCTACTCTCTCCGGCCAGATGCCCTTCTCCAGCTGCCCGCGTGTGGCAGCTCTACCCC  
 GGGAGGAAGTACCGAGCAGCGGATTCCAGTTTCTGGGAAATAGCAGACCATTGGCCAGATTGACTCTATG  
 GAGCAGCTGGTCCTAACATATCAGCCGGAGAGGAAAGAC**TGA**  
**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT  
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_199206
Insert Size:	387 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.


[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u>NM_199206.2</u>
<b>RefSeq Size:</b>	1738 bp
<b>RefSeq ORF:</b>	387 bp
<b>Locus ID:</b>	9623
<b>Cytogenetics:</b>	14q32.13
<b>MW:</b>	14.8 kDa
<b>Gene Summary:</b>	Enhances the phosphorylation and activation of AKT1 and AKT2.[UniProtKB/Swiss-Prot Function]