

Product datasheet for **SC307934**

SLC43A3 (NM_199329) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC43A3 (NM_199329) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLC43A3
Synonyms:	EEG1; ENBT1; FOAP-13; PRO1659; SEEEG-1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC307934 representing NM_199329.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGGGCCAGGGCCTGCCCTGCACGTGGCCACACTGCTGACTGGGCTGCTGGAATGCCTGGGCTTT
GCTGGCGTCCTCTTTGGCTGGCCTTCACTAGTGTCTTCAAGAATGAAGATTACTTTAAGGATCTG
TGTGGACCAGATGCTGGGCCGATTGGCAATGCCACAGGGCAGGCTGACTGCAAAGCCAGGATGAGAGG
TTCTCACTCATCTTCACCTGGGGTCTTCATGAACAACCTTCATGACATTCCTCCACTGGCTACATCTTT
GACCGGTTCAAGACCACCGTGGCAGCCTCATAGCCATATTTTTCTACACCACCGCCACACTCATCATA
GCCTTACCTCTGCAGGCTCAGCCGTGCTGCTCTTCTGGCCATGCCAATGCTCACCATTGGGGGAATC
CTGTTTCTCATCACCAACCTGCAGATTGGGAACCTATTTGGCCAACACCGTTCGACCATCATCACTCTG
TACAATGGAGCATTGACTCTTCTCGGCAGTCTTCTTATTATTAAGCTTCTTTATGAAAAGGCATC
AGCCTCAGGGCCTCTTCATCTTCTGTCTGCAGTACCTGGCATGTAGCACGCACTTCTCTCTG
ATGCCCCGGGGCACATCCATACCCACTGCCCCCAACTACAGCTATGGCCTGTGCCCTGGGAATGGC
ACCACAAAGGAAGAGAAGGAAACAGCTGAGCATGAAAACAGGGAGCTACAGTCAAAGGAGTTCCTTTCA
GCGAAGGAAGAGACCCAGGGGCAGGGCAGAAGCAGGAACCTCCGCTCCTTCTGGAGCTACGCTTCTCT
CGGCGCTTTGCTGGCACCTGGTGTGGCTGTCTGTGATACAGTTGTGGCACTACCTTTCATTGGCACT
CTCAACTCCTTGCTGACCAACATGGCCGGTGGGGACATGGCACGAGTCAACACCTACACAAATGCCTTT
GCCTTCACTCAGTTCGGAGTGTGTGTGCCCTTGGAAATGGCCTGCTCATGGACCGGCTTAACAGAAAG
TACCAGAAGGAAGCAAGAAAGACAGGTTCTCCACTTTGGCGTGGCCCTCTGCTCGACGGTGCCTTCG
CTGGCCCTGACATCCCTGCTGTGCCTGGGCTTCGCCCTCTGTGCCTCAGTCCCCATCCTCCCTCTCCAG
TACCTCACCTTCACTCTGCAAGTGATCAGCCGCTCCTTCTCTATGGGAGCAACGCGCCCTTCTCACC
CTTGCTTTCCCTCAGAGCACTTTGGCAAGCTTTGGGCTGGTATGGCCTTGTGGCTGTGGTGTCT
CTGCTCCAGTTCCCATTTACCCCTCATCAAAGGCTCCCTTCAAGATGACCCATTTTACGTAATGTG
ATGTTTCACTGCTTGCATTTCTGACATTTCTCCACCCCTTTCTGGTATATCGGGAATGCCGTAATTGG
AAAGAAAGTCCCTCTGCAATTGCATAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_199329

Insert Size: 1476 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_199329.2](#)

RefSeq Size: 2789 bp

RefSeq ORF: 1476 bp

Locus ID: 29015

UniProt ID: [Q8NBI5](#)

Cytogenetics: 11q12.1

Protein Families: Transmembrane

MW: 54.5 kDa

Gene Summary: Putative transporter.[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (3) differs in the 5' UTR, compared to variant 1. Variants 1, 2, 3 and 4 encode the same isoform 1.