

Product datasheet for **SC307950**

ILDR2 (NM_199351) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ILDR2 (NM_199351) Human Untagged Clone
Tag:	Tag Free
Symbol:	ILDR2
Synonyms:	C1orf32; dj782G3.1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC307950 representing NM_199351.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGATAGGGTCTTGCTGAGGTGGATTTCTCTCTTCTGGCTAACAGCCATGGTCAAGGCCCTTCAGGTC
ACAGTGGCCCGACAAGAAGAAGGTGGCCATGCTCTTCCAGCCACTGTGCTTCGCTGCCACTTCTCAACA
TCCTCCCATCAGCTGCAGTTGTGCAGTGGAAAGTTCAAGTCTACTGCCAGGATCGCATGGGAGAATCC
TTGGGCATGTCTCTACCCGGGCCAATCTCTCAGCAAGAGAAACCTGGAATGGGACCCCTACTTGAT
TGTTTGACAGCAGGAGGACTGTTTCAGTAGTAGCTTCAAACAGGGCTCGACTGTACCCTGGGAGAT
TTCTACAGGGGCAGAGAGATCACGATTGTTTCATGATGCAGATCTCAAATGGAAAGCTTATGTGGGA
GACAGCGGACTCTATTACTGTATTATCACACCCAGATGACCTGGAGGGGAAAAATGAGGACTCAGTG
GAACTGCTGGTGTGGCAGGACAGGGCTGCTTGTGATCTCTTGCCAGTTTGTGTGGAGATTATG
CCAGAGTGGTGTGGTGGCCTGGTGTCTTGGCGTCTTCTCTTCTCGTCTGGTGGGATCTGC
TGGTGCCAGTGTGCCCTCACAGCTGCTGCTATGTCCGCTGCCATGCTGCCAGATTCCTGCTGC
TGCCCTCAAGCCTTGATGAAGCAGGAAAGCAGCAAAGCCGGGTACCTCCCTGTCTCCGGTGTG
CCCGGCCCTTACTCCATCCCCTGTGCCCTTTGGGAGGAGCCCCCTCATCTGGCATGCTGATGGACAAG
CCGCATCCACCTCCCTTGGCACCAAGTGACTCCACTGGAGGAAGCCACAGTGTTCCGAAAGGTTACCGG
ATCCAGGCTGACAAAGAGAGAGACTCCATGAAGTCTGTACTATGTTGAGAAGGAGTGGCTCAGTTT
GATCCAGCCAGAAGGATGAGAGGAGATATAACAACACCATCTCAGAAGTCCAGCTCCCTACATGAGGAG
GACAGCAATTTCCGCCAGTCTTCCATCAGATGAGAAGCAAGCAGTCCCTGTGTCTGGGACTTGGAG
AGCAATCCTGACTATTGGTCAGGTGTATGGGAGGCAGCAGTGGGCAAGCCGGGCCCTCAGCCATG
GAGTATAACAAAGAGGATCGAGAGAGCTTCAAGCACAGCCAGCCGCTCCAAGTCGGAGATGCTGTG
CGGAAGAAGTTCGCCACGGGGTCCCGCGTTCATGGACGAGCTGGCGCCTTCGCTGACTCCTAC
GGCCAGCGGCCCGCGGGCAGACGCAACAGTACAGAGCGCGGGCGGGAGCCGCTTCGAGCGCTCG
GAGTCGCGGGCGCACAGCGGCTTACCAGGACTCCTTGGAGGAGTACTACGGTCAGCGCAGCCGC
AGCCGCGAGCCCCTGACCGATGCTGACCGCGGCTGGGCTTACGCCCCGCGCCGAGACCCGCGGAG
GACGCGCACCTGCCGCGGCTGGTGGCCGACGCCAGGCCACCGCAACAAATACGACCACTCGTACCTG
GGCAGCGCGGGAGCGCCAGGCGCGCCGAGGGCGCCAGCCGCGGTGGCAGCCTGGAGACGCCATCC
AAGCGGAGCGCGAGCTCGGCCCGCGCAGCGCCTCTACTACGCTTGGTCGCGCCCGGCACCTACAAG
GCCGGCTCGTCGAGGACGACCAGGAGGACGCGTCCGACGACGCGCTGCCGCCCTACAGCGAGCTGGAG
CTGACCCGCGGCCGTCCTACCGCGGCCGACCTGCCCTACCACAGCAACTCGGAGAAGAAGAGGAAA
AAGGAGCCCGCAAGAAAACCAATGACTTTCAACCAGGATGTCCCTTGTGGTCTGA
AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
  
```

Restriction Sites: SgfI-RsrII

Plasmid Map: □

ACCN: NM_199351

Insert Size: 1920 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_199351.2</u>
RefSeq Size:	8127 bp
RefSeq ORF:	1920 bp
Locus ID:	387597
UniProt ID:	<u>Q71H61</u>
Cytogenetics:	1q24.1
Protein Families:	Transmembrane
MW:	71.2 kDa
Gene Summary:	May be involved in lipid homeostasis and ER stress pathways.[UniProtKB/Swiss-Prot Function]